

BATRIK Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long, ergonomic designed handle permits ease-of-use
- Extra long length of brush head (60 [\pm 5] mm in Single-Use and 75 [\pm 5] mm in Reusable)
- The brushes fibers are transparent, allowing for quick verification of its cleanliness
- The plastic composition used to fabricate the handle will not release any toxic components



BATRIK manufactures four types of flat brushes:

- **EXTRA RIGID** instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.
- **RIGID** the rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.
- **FLEXIBLE** This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics or hard to reach spots which are easier to clean using a flexible bristle.
- **EXTRA FLEXIBLE** The very long bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments (optics, etc.) as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach with shorter bristles, such as the accessory channel or tip deflection of pliable endoscopes



Extra Rigid Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40917	Reusable Flat brush (extra rigid)	220mm	75L x 20W x 10H

Packaging: Bag of 3



Rigid Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40918	Reusable Flat brush (rigid)	220mm	75L x 20W x 15H

Packaging: Bag of 3



Flexible Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40919	Reusable Flat brush (flexible)	220mm	75L x 20W x 15H

Packaging: Bag of 3



Extra Flexible Reusable / Autoclavable Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40920	Reusable Flat brush (extra flexible)	220mm	75L x 20W x 25H

Packaging: Bag of 3



Rigid and Flexible Single Use Flat Brush

Item #	Description	Length – mm	Brush Size – mm
50-40921	Single use flat brush (rigid)	220mm	60L x 20W x 15H
50-40922	Single use flat brush (flexible)	220mm	60L x 20W x 15H

Packaging: Bag of 20

POLYAMIDE FIBERS

