


BATRIK STERILIZATION Reusable Orthopedic Rasps Cleaning Brush

NEW!
Antibacterial

This reusable brush is used to clean orthopedic rasps and other similar instruments in orthopaedic surgery.

Features

- Antibacterial polyamide bristles
- Ergonomic steamable thermoplastic handle
- Reusable
- Long lasting and can withstand many sterilization cycles
- Latex free



Orthopedic Rasps Cleaning Brush

Item #	Total Length Inches / cm	Handle Length Inches / cm	Brush Size Inches / cm	Packaging
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	3 Units / Bag

ANTIBACTERIAL

The bristles are made of polyamide fibers which prevent or slow bacterial growth between two cleanings or brush sterilizations.