

LEVEL 2 PROTECTIVE GOWNS

LEVEL 2
NON-STERILE
ULTRASONICALLY WELDED

Our **Level 2** Disposable Protective Water Repellant Gowns are designed to protect health professionals and patients from cross contamination.

- Compliant to EN ISO 22610
- Compliant to EN ISO 20811 LIQUID BARRIER PERFORMANCE
- SMS 40 gr/m² soft non-woven spunbond fabric provides **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Non-Sterile



SURGMED EST. 1992.

All our Medical Devices are product liability insured

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-202-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-202-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-202-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-202-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"

PACKAGING: 50 GOWNS / BOX | FOR SALE IN USA & CANADA



CONTACT US:

BATRIK MEDICAL MANUFACTURING INC.
850 Halpern Ave., (Dorval) Montreal QC H9P 1G6 CANADA
T.F.: 1-800-661-5432 | T.: +1 514-631-7988 | F.: + 1-514-631-9083 | info@surgmed.com | www.batrik.com