

# LEVEL 4 PROTECTIVE LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

**LEVEL 4**  
**NON-STERILE**  
ULTRASONICALLY WELDED  
& REINFORCED SEAMS

Our **Level 4** Disposable Protective Laminated Gowns provide high performance for use in high risk situations and/or when infectious diseases are suspected.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance** and blood and viral penetration resistance
- 56 gr/m<sup>2</sup> soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Non-Sterile**



INSIDE VIEW OF  
ULTRASONICALLY WELDED  
& REINFORCED SEAMS

**SURGMED EST. 1992.**

All our Medical Devices are product liability insured

## PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
<b>SF-402-S</b>	<b>S</b>	110 cm / 43"	55 cm / 22"	140 cm / 55"
<b>SF-402-M</b>	<b>M</b>	115 cm / 45"	58 cm / 23"	145 cm / 57"
<b>SF-402-L</b>	<b>L</b>	127 cm / 50"	60 cm / 24"	150 cm / 59"
<b>SF-402-XL</b>	<b>XL</b>	135 cm / 53"	62 cm / 24"	150 cm / 59"
<b>SF-402-XXL</b>	<b>XXL</b>	140 cm / 55"	64 cm / 25"	155 cm / 61"
<b>SF-402-XXXL</b>	<b>XXXL</b>	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX | FOR SALE IN CANADA



### CONTACT US:

BATRIK MEDICAL MANUFACTURING INC.  
850 Halpern Ave., (Dorval) Montreal QC H9P 1G6 CANADA  
T.F.: 1-800-661-5432 | T.: +1 514-631-7988 | F.: +1-514-631-9083 | info@surgmed.com | www.batrik.com