

CHEMOTHERAPY PROTECTIVE LAMINATED ULTRASONICALLY WELDED GOWNS

LEVEL 3
NON-STERILE
ULTRASONICALLY WELDED

Our **Level 3** Disposable Chemotherapy Protective Laminated Gowns are impermeable to **chemotherapy drugs** and other fluids, offering outstanding protection.

- Fabric Compliant to EN 20811 and EN ISO 22610
- Tested and compliant to AAMI PB70 - Level 3 Liquid Barrier Performance
- Tested and compliant to ASTM F739-12
- Chemo Resistance to **12 CHEMOTHERAPY DRUGS: BICNU, Adriamycin, Cytosan, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin**
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- Ultrasonically welded seams
- Non-Sterile



SURGMED EST. 1992.
All our Medical Devices are product liability insured

PRODUCT SELECTION

*ITEM #	**ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-502-S	SF-602-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-502-M	SF-602-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-502-L	SF-602-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-502-XL	SF-602-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-502-XXL	SF-602-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-502-XXXL	SF-602-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

* PACKAGING SF-502: 40 INDIVIDUALLY BAGGED GOWNS / BOX | FOR SALE IN CANADA

** PACKAGING SF-602: 50 BULK PACKED GOWNS / BOX | FOR SALE IN CANADA



CONTACT US:

BATRIK MEDICAL MANUFACTURING INC.
850 Halpern Ave., (Dorval) Montreal QC H9P 1G6 CANADA
T.F.: 1-800-661-5432 | T.: +1 514-631-7988 | F.: + 1-514-631-9083 | info@surgmed.com | www.batrik.com