

# STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED GOWNS

**LEVEL 3  
STERILE**  
ULTRASONICALLY WELDED

Our **Level 3 Disposable Sterile Surgical Laminated Gowns** provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 3 Liquid Barrier Performance**
- Tested and compliant to **ASTM F739-12**
- Chemo Resistance to **12 CHEMOTHERAPY DRUGS: BiCNU, Adriamycin, Cytosin, Taxol, Toposar, 5-Flu, Thiotepa, Cisplatin, DTIC, Mitomycin C, Oncovin, and Amethopterin**
- 56 gr/m<sup>2</sup> fabric for maximum **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded seams**
- **Sterile and aseptically packed**

Individually packed, each gown comes in a set designed to be used during surgery to **reduce contamination and infection** between the patient and medical team.

#### Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 3 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")



**SURGMED EST. 1992.**

All our Medical Devices are product liability insured

## PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
<b>41.001.01</b>	<b>S</b>	110 cm / 43"	55 cm / 22"	140 cm / 55"
<b>41.002.01</b>	<b>M</b>	115 cm / 45"	58 cm / 23"	145 cm / 57"
<b>41.003.01</b>	<b>L</b>	127 cm / 50"	60 cm / 24"	150 cm / 59"
<b>41.004.01</b>	<b>XL</b>	135 cm / 53"	62 cm / 24"	150 cm / 59"
<b>41.005.01</b>	<b>XXL</b>	140 cm / 55"	64 cm / 25"	155 cm / 61"
<b>41.006.01</b>	<b>XXXL</b>	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 40 individually bagged sets / box | FOR SALE IN CANADA



#### CONTACT US:

Advanced Surgi-Pharm Inc.  
850 Halpern Ave., (Dorval) Montreal QC H9P 1G6  
TF: 1-800-661-5432 | T.: +1-514-631-7988 | F.: +1-514-631-9083 | E.: info@surgmed.com | www.surgmed.com