

STERILE SURGICAL LAMINATED ULTRASONICALLY WELDED & REINFORCED GOWNS

**LEVEL 4
STERILE**
ULTRASONICALLY WELDED
& REINFORCED SEAMS

Our **Level 4 Disposable Sterile Surgical Laminated Gowns** provide high performance with a fluid resistant barrier protection.

- Fabric Compliant to **EN 20811** and **EN ISO 22610**
- Tested and compliant to **AAMI PB70 - Level 4 Liquid Barrier Performance** and blood and viral penetration resistance
- 56 gr/m² soft non-woven spunbond fabric coated with a polyethylene film provides maximum **protection, breathability, and comfort**
- Tested and compliant to **ASTM F739-12**
- Knitted soft wrist cuffs
- Tie closure
- **Ultrasonically welded & reinforced seams**
- **Sterile and aseptically packed**

Individually packed, each gown comes in a set designed to be used during surgery to **reduce contamination and infection** between the patient and medical team.

Each set includes:

- 1 x Sterile Surgical Laminated Gown - Level 4 (Size S - XXXL)
- 1 x Blue Sterilization Wrap (60 cm x 60 cm / 16" x 16")
- 1 x White Cellulose Hand Towel (30 cm x 40 cm / 12" x 16")



INSIDE VIEW OF
ULTRASONICALLY WELDED
& REINFORCED SEAMS

SURGMED EST. 1992.

All our Medical Devices are product liability insured

PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
46.001.01	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
46.002.01	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
46.003.01	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
46.004.01	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
46.005.01	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
46.006.01	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 40 INDIVIDUALLY BAGGED SETS / BOX | FOR SALE IN CANADA



CONTACT US:

Advanced Surgi-Pharm Inc.
850 Halpern Ave., (Dorval) Montreal QC H9P 1G6
TF: 1-800-661-5432 | T.: +1-514-631-7988 | F.: +1-514-631-9083 | E.: info@surgmed.com | www.surgmed.com