


Your First ™

## NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS



Presenting our enhanced Non-Enzymatic Bedside Sponge Kits. Each kit includes sponges pre-treated with a free-rinsing, neutral pH detergent blend from RightSwipe™, ensuring an effective pre-cleaning and surface cleaning solution for medical instruments. This kit stands as Your First ™ step in the medical instrument cleaning process, showcasing the advanced benefits for superior instrument care. Our enhanced bedside kits will provide you with a streamlined workflow.



FLAT SPONGE



TUBULAR SPONGE



FREESTANDING POUCH



ENVIRO-SACHET™  
NON-ENZYMATIC CLEANER



DISPENSER BOX

## PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING	CASE DIMENSIONS
717-122	PRE-TREATED TUBULAR SPONGE Add 100-500 ml of water to activate non-enzymatic detergent. Comes in plastic bag.1	15 / Box (4 Boxes / Case)	21" x 14" x 13.5"
717-222	PRE-TREATED FLAT SPONGE Add 100-500 ml of water to activate non-enzymatic detergent. Comes in plastic bag.1	25 / Box (4 Boxes / Case)	21" x 14" x 13.5"
717-133	READY TO USE TUBULAR SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
717-233	READY TO USE FLAT SPONGE Comes pre-treated in 500 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
717-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" x 14" x 13.5"
717-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes pre-treated in 100 ml of non-enzymatic solution in freestanding pouch.	12 / Box (4 Boxes / Case)	21" x 14" x 13.5"
717-155	READY TO USE TUBULAR SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
717-255	READY TO USE FLAT SPONGE Comes pre-treated in 220 ml of non-enzymatic solution in freestanding pouch.	24 / Case	14" x 12" x 9"
717-277	NON-ENZYMATIC CLEANER Comes with 1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit.	1 Kit	1 Kit

Our stainless steel bowl can be sold separately and used as a basin for easy and convenient sponge activation.1