


# BATRIK VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES  
AND 4 COLORS



Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.


Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

## PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" 0.2 cm x 40 cm	Mini-Red
10-1001	0.075" X 16" 0.2 cm x 40 cm	Mini-Blue
10-1002	0.075" X 16" 0.2 cm x 40 cm	Mini-Yellow
10-1003	0.075" X 16" 0.2 cm x 40 cm	Mini-White
10-1010	0.120" x 16" 0.3 cm x 40 cm	Maxi-Red
10-1011	0.120" x 16" 0.3 cm x 40 cm	Maxi-Blue
10-1012	0.120" x 16" 0.3 cm x 40 cm	Maxi-Yellow
10-1013	0.120" x 16" 0.3 cm x 40 cm	Maxi-White
10-1101	0.040" x 16" 0.1 cm x 40 cm	Extra-Mini-Blue
10-1111	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-Blue
10-1113	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-White



- Manufactured of medical grade silicone



- Do not twist on themselves when unpacked



- Easy to separate, eliminating the need the need for cardboard inserts



- Sterile, radiopaque, and latex free

PACKAGING: 2 ties per sterile pouch – 10 pouches per box