


BATRIK ENDOSCOPY BO2™ BITE BLOCKS WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

BO2
(FRONT SIDE)

BO2
(BACK SIDE)



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
BO2	Single-Use	18mm / 54FR, includes strap	100 / Box
BO2-N	Single-Use	18mm / 54FR, without strap	100 / Box

BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

- Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
B1	Single-Use Pediatric	16 mm / 48FR, includes strap	100 / Box
B1-25	Single-Use Pediatric	16 mm / 48FR, includes strap	25 / Box
B1-N	Single-Use Pediatric	16 mm / 48FR, without strap	100 / Box
B3	Single-Use	20 mm / 60FR, includes strap	100 / Box
B3-N	Single-Use	20 mm / 60FR, without strap	100 / Box
B3-SA	Single-Use	20 mm / 60FR, with strap <u>attached</u> to one side	100 / Box
B3-NDR	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box
B3-NDR-N	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box
B4	Single-Use	19 mm / 60FR	100 / Box
B4T	Single-Use	19 mm / 60FR, includes strap	100 / Box
B4T-N	Single-Use	19 mm / 60FR, without strap	100 / Box
B6	Reusable	19 mm / 60FR	100 / Box
B6T	Reusable	19 mm / 60FR, includes strap	100 / Box
B6T-N	Reusable	19 mm / 60FR, without strap	100 / Box
B6-25	Reusable	19 mm / 60FR	25 / Box
B6T-25	Reusable	19 mm / 60FR, includes strap	25 / Box
B6T-N-25	Reusable	19 mm / 60FR, without strap	25 / Box
BS	Single-Use	Bite-Block straps 21.7" x .79" / 55 cm x 2 cm	100 / Box



B5



B1



B3



B4



B6