


BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




NEW!

18-CP06

COLOUR	BLUE	RED	PINK
USAGE	SINGLE-USE	REUSABLE <small>withstands repeated sterilization up to 134°C</small>	REUSABLE <small>withstands repeated sterilization up to 134°C</small>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors feature foam that is 0.25" thicker than competitor products. They are especially useful with heavier instrument sets.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
18-CP02 **	Corner Pockets™ Single-Use, Blue	240 / Bag
18-EVAL2 **	Corner Pockets™ Single-Use, Blue	24 / Bag
18-CP07	Corner Pockets™ Reusable, Red	80 / Bag
18-CORNER	Corner Pockets™ Reusable, Red	320 / Box (4 Bags of 80 pcs)
18-EVAL7	Corner Pockets™ Reusable, Red	24 / Bag
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.