




# BATRIK INSTRUMENT PROTECTORS

Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization
- Multiple color-coded sizes
- Latex Free
- Non-Sterile



## PRODUCT SELECTION

### REGULAR & PERFORATED / CLEAR



RC-07



PC-07

ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
<b>RC-01</b>	<b>PC-01</b>	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
<b>RC-02</b>	<b>PC-02</b>	Clear	0.075" x 0.75" (2 mm x 19 mm)
<b>RC-03</b>	<b>PC-03</b>	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
<b>RC-04</b>	<b>PC-04</b>	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
<b>RC-05</b>	<b>PC-05</b>	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
<b>RC-06</b>	<b>PC-06</b>	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
<b>RC-07</b>	<b>PC-07</b>	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
<b>RC-08</b>	<b>PC-08</b>	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
<b>RC-09</b>	<b>PC-09</b>	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
<b>RC-10</b>	<b>PC-10</b>	Clear	Assorted

PACKAGING: Bag of 100

### REGULAR & PERFORATED / TINTED



RT-07



PT-07

<b>RT-01</b>	<b>PT-01</b>	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
<b>RT-02</b>	<b>PT-02</b>	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
<b>RT-03</b>	<b>PT-03</b>	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
<b>RT-04</b>	<b>PT-04</b>	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
<b>RT-05</b>	<b>PT-05</b>	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
<b>RT-06</b>	<b>PT-06</b>	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
<b>RT-07</b>	<b>PT-07</b>	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
<b>RT-08</b>	<b>PT-08</b>	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
<b>RT-09</b>	<b>PT-09</b>	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
<b>RT-10</b>	<b>PT-10</b>	Tinted Assorted	Assorted

PACKAGING: Bag of 100

### REGULAR & PERFORATED / SOLID



RS-07



PS-07


<b>RS-01</b>	<b>PS-01</b>	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
<b>RS-02</b>	<b>PS-02</b>	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
<b>RS-03</b>	<b>PS-03</b>	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
<b>RS-04</b>	<b>PS-04</b>	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
<b>RS-05</b>	<b>PS-05</b>	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
<b>RS-06</b>	<b>PS-06</b>	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
<b>RS-07</b>	<b>PS-07</b>	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
<b>RS-08</b>	<b>PS-08</b>	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
<b>RS-09</b>	<b>PS-09</b>	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
<b>RS-10</b>	<b>PS-10</b>	Solid Assorted	Assorted

PACKAGING: Bag of 100

ITEM #	DIMENSIONS
<b>OP-11</b>	0.25" x 0.75" (6.35 mm x 19.05 mm)
<b>OP-12</b>	0.375" x 1" (9.53 mm x 25.4 mm)
<b>OP-13</b>	0.5" x 1" (12.7 mm x 25.4 mm)
<b>OP-14</b>	0.75" x 1" (19.05 mm x 25.4 mm)
<b>OP-15</b>	1" x 1" (25.4 mm x 25.4 mm)
<b>OP-16</b>	1.25" x 1.25" (31.75 mm x 31.75 mm)
<b>OP-17</b>	1.5" x 1.5" (38.1 mm x 38.1 mm)
<b>OP-18</b>	1.75" x 1.25" (44.45 mm x 31.75 mm)
<b>OP-20</b>	Assorted

**PACKAGING:** Bag of 50

**OSTEOTOME PROTECTORS**



ITEM #	DIMENSIONS	PACKAGING
<b>RGS-01</b>	1.5" (38 mm)	55 / Bag
<b>RGS-02</b>	2" (51 mm)	45 / Bag
<b>RGS-03</b>	3" (76.2 mm)	40 / Bag


**SINGLE-TIP RETRACTOR PROTECTORS**



ITEM #	COLOR	SIZE
<b>RGM-01</b>	Pink	Small
<b>RGM-02</b>	Clear	Small
<b>RGM-03</b>	Green	Medium
<b>RGM-04</b>	Clear	Medium

**PACKAGING:** Bag of 50

**WEITLANER RETRACTOR PROTECTORS**



ITEM #	COLOR	DIMENSIONS
<b>DG-01</b>	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
<b>DG-02</b>	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
<b>DG-03</b>	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
<b>DG-04</b>	Red	0.125" x 1" (3.18 mm x 25.4 mm)
<b>DG-05</b>	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
<b>DG-07</b>	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
<b>DG-10</b>	Assorted	Assorted

**PACKAGING:** Bag of 50

**PERFORATED DOUBLE PROTECTORS**

