




# BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



**FRONT SIDE**  
Pins keep instruments securely in place



**BACK SIDE**  
Numerous holes and elevated nubs provide superior venting and reduce water collection

## PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
<b>13-2001</b>	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	1 unit
<b>13-2005</b>	13" x 18" (33 cm x 45.72 cm) 11.97 mm high from base of mat to end of pin 15.53 mm high from base of elevated nub to end of pin	5 units

GOOD TO KNOW:  
CUSTOM SIZES  
ARE AVAILABLE  
UPON REQUEST