

Your First ™

PrAmLi Power³ Enzymes™

ENZYMATIC DETERGENTS

All **PrAmLi Power³ Enzymes™** detergents are multi-enzyme pre-soak cleaners for medical equipment and instrumentation. Our proprietary unique blend of enzymes is effective in the removal of all proteins, blood, fatty tissue, and more. It will also remove stains from surgical scrubs, linens, and medical clothing.



Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



PRODUCT SELECTION

ITEM #	DILUTION DETERGENT / WATER	PACKAGING
LOW FOAM DETERGENT		
SS-21-1001	4 to 8 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-1004	4 to 8 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
MILD FOAM DETERGENT		
SS-21-2001	4 to 8 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-2004	4 to 8 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
SUPER CONCENTRATE		
SS-21-3001	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-3004	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	4 x 3.8 L (4 Gallons)

NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.

GOOD TO KNOW:

OUR LOW FOAM DETERGENT MAY BE USED IN WASHERS/ ULTRASONIC CLEANERS. WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER/ULTRASONIC CLEANER AS MACHINES/ SETTINGS VARY.