

BATRIK

# ARC™ WIRE STYLE BRUSHES

## SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

### BRUSH SELECTION TIPS:

#### INSTRUMENT SPECIFICATIONS

##### What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

##### What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that

are too large will bend too much and not clean effectively.

#### BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.

In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At Batrik, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

#### CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

## RIGID



## FLEXIBLE



**12" / 30 CM SINGLE-END WIRE BRUSH**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-31502-S</b>	<b>41-31502</b>	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-32503-S</b>	<b>41-32503</b>	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-32504-S</b>	<b>411-32504</b>	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-32505-S</b>	<b>411-32505</b>	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-35005-S</b>	<b>41-35005</b>	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35006-S</b>	<b>41-35006</b>	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-35008-S</b>	<b>41-35008</b>	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35010-S</b>	<b>41-35010</b>	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35012-S</b>	<b>41-35012</b>	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-35015-S</b>	<b>41-35015</b>	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-35020-S</b>	<b>41-35020</b>	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

**18" / 45 CM SINGLE-END WIRE BRUSH**

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-46502-S</b>	<b>41-46502</b>	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-47503-S</b>	<b>41-47503</b>	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-47504-S</b>	<b>411-47504</b>	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-47505-S</b>	<b>411-47505</b>	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-50005-S</b>	<b>41-50005</b>	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50006-S</b>	<b>41-50006</b>	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-50008-S</b>	<b>41-50008</b>	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50010-S</b>	<b>41-50010</b>	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50012-S</b>	<b>41-50012</b>	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-50015-S</b>	<b>41-50015</b>	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-50020-S</b>	<b>41-50020</b>	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

**24" / 60 CM SINGLE-END WIRE BRUSH**



	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-61502-S</b>	<b>41-61502</b>	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	<b>41-62503-S</b>	<b>41-62503</b>	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	<b>411-62504-S</b>	<b>411-62504</b>	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	<b>411-62505-S</b>	<b>411-62505</b>	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-65005-S</b>	<b>41-65005</b>	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65006-S</b>	<b>41-65006</b>	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	<b>41-65008-S</b>	<b>41-65008</b>	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65010-S</b>	<b>41-65010</b>	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65012-S</b>	<b>41-65012</b>	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-65015-S</b>	<b>41-65015</b>	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-65020-S</b>	<b>41-65020</b>	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

**30" / 75 CM SINGLE-END WIRE BRUSH**


	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
<b>FLEXIBLE</b> BRUSH & SHEATH	<b>41-76502-S</b>	<b>41-76502</b>	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	<b>41-77503-S</b>	<b>41-77503</b>	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	<b>411-77504-S</b>	<b>411-77504</b>	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	<b>411-77505-S</b>	<b>411-77505</b>	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
<b>RIGID</b> BRUSH & SHEATH	<b>41-80005-S</b>	<b>41-80005</b>	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80006-S</b>	<b>41-80006</b>	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	<b>41-80008-S</b>	<b>41-80008</b>	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80010-S</b>	<b>41-80010</b>	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80012-S</b>	<b>41-80012</b>	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	<b>41-80015-S</b>	<b>41-80015</b>	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	<b>41-80020-S</b>	<b>41-80020</b>	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



## FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



## RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

### BRUSH SELECTION SYSTEM (Actual Size)

2mm 0.078"	3 mm 0.118"	4 mm 0.157"	5 mm 0.196"	6 mm 0.236"	8 mm 0.314"	10 mm 0.394"	12 mm 0.472"	15 mm 0.591"	20 mm 0.787"