



BATRIK SAVEFACE™ FACE SHIELDS

Our Batrik SaveFace™ Face Shields (Visors) are lightweight and allow for good peripheral vision while acting as a physical barrier against splashes. Can be used for extended periods due to comfort-focused design and superior quality materials.

Choose amongst non-draped and draped models:

- Our classic **non-draped** models close the gap between the shield and the chest bone when head is tilted forward for added protection.
- Our **draped** models offer full-length coverage extending from the top of the eyebrow to below the chest, covering the anterior neck area as per CDC recommendations.



- Self-shaping foam headpieces mold and cushion forehead for a comfortable and secure fit (per CDC recommendations)
- Shielded air vents maximize airflow for enhanced comfort and prevent fogging
- Excellent Distortion-Free Optical Clarity with flat viewing area reduces distortion and glare
- Features anti-fog treatment for clear view
- Wide visors provide peripheral protection to help prevent splashes from going into the eyes
- High Quality PET
 - SF-100/SF-180: 7.5 mil (0.19 mm) thick
 - SF-240/SF-8240: 9.84 mil (0.25 mm) thick
- Our patent-pending drape models come equipped with integrated ANSI/AAMI PB70 Level 3 compliant splash drape barriers for added protection

ASSEMBLED IN CANADA

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-100	SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	100 / Box
SF-180	SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	100 / Box
SF-240	BATRIK SAVEFACE™ no drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm	9" (23 cm)	240 / Case (8 Boxes of 30 Pcs)
SF-8240	BATRIK SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Shield Width: 13" / 33 cm Drape Length: 6" / 15 cm	15" (38 cm)	240 / Case (8 Boxes of 30 Pcs)