

# BATRIK REUSABLE INSTRUMENT CLEANING BRUSHES


Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

### REUSABLE JAR/BOTTLE BRUSH



### STERILIZER DRAIN BRUSH



Blue Tip  
for added safety

Heavy Duty  
Stainless Wire

Ergonomic  
Handle

## PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
<b>50-10295</b>	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
<b>50-10224</b>	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm

PACKAGING: Bag of 3