

BATRIK HYSTEROSONOGRAPHY CATHETER

Hystero-sonography (also known as sonohysterography) can now be performed safely, simply and atraumatically with the gentle infusion of saline into the uterus with the use of our H/S Catheter.


- Available in two sizes: 5FR and 7FR
- Minimally invasive and ideal for office use
- Gentle placement sheath included (decreasing the need for a tenaculum) to minimize patient discomfort
- Flexible and easy-to-use
- Total length of 410mm with a working length of 280 mm
- Latex Free

GOOD TO KNOW: ACHIEVE
IMPROVED IMAGERY
AND DIAGNOSIS WITH
SONOHYSTEROGRAMS

“ We depend upon the reliability, precision and functionality of the Batrik H/S Catheter in our clinical and research operations.

Having performed more than 100,000 sonohysterography procedures, we can unequivocally say that this is the best HSG catheter that we have ever used

”
-Dr. Alex Hartman



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
51-5005	5 FR H/S Catheter	10 / box
51-5007	7 FR H/S Catheter	10 / box