



COMPLETE PRODUCT CATALOGUE

Established in 1997. Proudly Canadian. 🍁
Sold in 80 countries worldwide.



ABOUT US



Proudly Canadian and established in 1997, we specialize in the manufacturing of enzymatic solutions offered in different formats as well as the fabrication of silicone and plastic products. Our products are sold in 80 countries worldwide and used primarily in surgery, endoscopy, and central sterilization areas.

Our in-house controlled environment production facility in Montreal, Canada adheres to strictly enforced quality and regulatory control systems. We are registered as a device manufacturing establishment with the United States FDA and we are in full compliance with CE and Health Canada Regulations. We continually invest in upgrading our plant, machinery, and processes to ensure high quality products at high volumes. Most products are available off-the-shelf.

At Batrik, we are committed to innovation and providing our customers with the finest quality products backed by world-class customer care.

SURGICAL PRODUCTS & GENERAL ACCESSORIES









4	SUTURE CLAMP TAGS
5	CLAMP JAW COVERS
6	VASCULAR SILICONE TIES
7	CAUTERY TIP CLEANER
8	SKINS™ Surgical MULTI-TIP MarkerS
9	PIN GUARDS
10	SAFE-SOAK™ SPAC & MATS
11	MAGNETIC MAT
12	TOURNY™ + TOURNIQUETS
14	SAVEFACE™ SHOE COVERS
15	WTFOG™ Reusable Anti-Fog Cloths
17	STACKABLE STEP STOOLS
18	BLUE CLOUD-GEL™ ANTI-FATIGUE MATS
19	BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

ENDOSCOPY PRODUCTS

20	G.O.L.F.F.™ ANTI-FOG SOLUTION
21	BITE BLOCKS Endoscopy Bo2™ / Bitestrong™
23	COUVRETTE™ VALVE COVERS
24	DAYBRUSH™ CHANNEL BRUSHES
25	ARC™ CHANNEL BRUSHES
28	VALVE / CONTROL BRUSHES
29	V-BRUSHES

INFECTION CONTROL & STERILE PROCESSING

30	SILICONE PIN MAT
31	CORNER POCKETS™ & SOFT POCKETS™
32	INSTRUMENT PROTECTORS
35	PrAmLi Power³ LUBE™ INSTRUMENT LUBRICANT
36	S.T.I.C.S.
37	ID TAPE
38	TAPE-B-GONE™ TAPE REMOVER
39	PADLOCKS SECURITY TAGS
43	SAVEFACE™ FACE SHIELD
44	SHIELD™ FACE SHIELD
45	SCRUB PAD

46	Your First  ENZYMATIC DETERGENTS
47	Your First  ENVIRO-SACHET™
48	Your First  FOAMING SPRAY
49	Your First  F1RST WIPES
50	Your First  DRY SPONGES
51	Your First  ENZYMATIC BEDSIDE SPONGE KITS
52	Your First  NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS
53	Your First  TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

BRUSHES

54	ARC™ WIRE STYLE BRUSHES
58	DOUBLE ARC™ WIRE BRUSHES Asymmetrics / Symmetric
61	IN-BRUSH™ TWISTED WIRE BRUSHES
63	STAINLESS STEEL ARC™ MICRO BRUSHES
64	INSTRUMENT CLEANING FLAT BRUSHES
65	ORTHOPEDIC RASPS CLEANING BRUSH
66	PRECISION BRUSH
67	JAR / BOTTLE & DRAIN BRUSHES
68	ACETABULAR REAMER BRUSH
69	ALL PURPOSE REUSABLE BRUSHES
70	AUTOCLAVE CLEANING BRUSH KITS

WOMEN'S HEALTH


71	HYSTEROSONOGRAPHY CATHETER
----	----------------------------

BATRIK SUTURE CLAMP TAGS


Non-toxic, radiopaque, and latex free, our Suture Clamp Tags are designed to protect and safely hold sutures without any damage.

Additionally, its foam base, affixed to a sterile area with the clamps, allows for easy insertion.

- Packed in a convenient foam base with easy peel-off adhesive backing, even while wearing gloves
- Available in sterile and non-sterile
- Available in 2, 3, or 5 pairs
- Pre-counted



OUR "EASY PEEL CORNER" IS A FEATURE UNIQUE TO BATRIK



PRODUCT SELECTION

ITEM #	STERILE	ITEM#	NON-STERILE	COLOR	DIMENSIONS	PACKAGING
11-1002		911-1002		Yellow	Standard ø (0.06" x 0.5" / 0.15 cm x 1.27 cm)	2 pairs / foam base
11-1003		911-1003		Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	3 pairs / foam base
11-1005		911-1005		Yellow	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1025		911-1025		Assorted	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1045		911-1045		White	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1065		911-1065		Blue	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-1085		911-1085		Red	Standard ø (0.06" x 0.5" L / 0.15 cm x 1.27 cm)	5 pairs / foam base
11-3002		911-3002		Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	2 pairs / foam base
11-3003		911-3003		Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	3 pairs / foam base
11-3005		911-3005		Pink	Large ø (0.09" x 0.75" L / 0.23 cm x 1.91 cm)	5 pairs / foam base

STERILE PACKAGING: 1 Foam Base / Pouch – 5 Pouches / Box


NON-STERILE PACKAGING: Sold per foam base

BATRIK CLAMP JAW COVERS

Manufactured from medical grade silicone, our Clamp Jaw Covers reduce trauma to vessels requiring occlusion, particularly when a hemostatic clamp is utilized. The thread-like indents allow for gentle, constant and non-slip tension.



GOOD TO KNOW: OUR CLAMP JAW COVERS CONVERT ALMOST ANY INSTRUMENT TO ATRAUMATIC



- Gamma sterilized
- Available in 3 colors and 3 widths to fit most clamp sizes
- Radiopaque and latex free

PRODUCT SELECTION

ITEM # STERILE	ITEM # NON-STERILE	DIMENSIONS	COLOR
12-1000	912-1000	0.1" X 5" / 0.25 cm X 12.7 cm	Red
12-1001	912-1001	0.2" X 5" / 0.51 cm X 12.7 cm	White
12-1002	912-1002	0.3" X 5" / 0.76 cm X 12.7 cm	Blue

STERILE PACKAGING: 2 clamps per sterile pouch - 10 pouches per box

NON-STERILE PACKAGING: 20 clamps per bag (also available 100 units / bag)

BATRIK VASCULAR SILICONE TIES

AVAILABLE IN 4 SIZES
AND 4 COLORS



Our Vascular Silicone Ties permit safe and smooth retraction, occlusion and identification of vessels during surgical procedures.

Unlike clamps and retractors which can pinch the vessels flat, or even crack them, the BATRIK elliptical shape allows for the gentle handling of arteries, veins, nerves, tendons, spermatic cords, ureters, or whenever delicate retraction / occlusion is required.

PRODUCT SELECTION

ITEM #	DIMENSIONS	COLOR
10-1000	0.075" X 16" 0.2 cm x 40 cm	Mini-Red
10-1001	0.075" X 16" 0.2 cm x 40 cm	Mini-Blue
10-1002	0.075" X 16" 0.2 cm x 40 cm	Mini-Yellow
10-1003	0.075" X 16" 0.2 cm x 40 cm	Mini-White
10-1010	0.120" x 16" 0.3 cm x 40 cm	Maxi-Red
10-1011	0.120" x 16" 0.3 cm x 40 cm	Maxi-Blue
10-1012	0.120" x 16" 0.3 cm x 40 cm	Maxi-Yellow
10-1013	0.120" x 16" 0.3 cm x 40 cm	Maxi-White
10-1101	0.040" x 16" 0.1 cm x 40 cm	Extra-Mini-Blue
10-1111	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-Blue
10-1113	0.200" x 22" 0.5 cm x 56 cm	Extra-Maxi-White



- Manufactured of medical grade silicone



- Do not twist on themselves when unpacked



- Easy to separate, eliminating the need the need for cardboard inserts



- Sterile, radiopaque, and latex free

PACKAGING: 2 ties per sterile pouch – 10 pouches per box

BATRIK CAUTERY TIP CLEANER

Our Cautery Tip Cleaner assists in removing material from cautery tips for both mono-polar and bi-polar electrodes.

- Foam pad with integrated abrasive material
- Adheres to disposable and non-disposable drapes
- Radiopaque
- Available exclusively in Non-Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-10000	Cautery Tip Cleaner Non-Sterile	1.8" x 1.8" / 4.5 cm x 4.5 cm	2,200 /Box

BATRIK SKINS™ SURGICAL MULTI-TIP SKIN MARKERS

Our SKINS™ Surgical Multi-Tip Skin Markers are designed with a beveled tip to enable the creation of fine or bold tracing lines by holding the pen at different angles.

- Available with or without 6 labels and a flexible 6" calibrated (in/mm) ruler
- Available in non-sterile for pre-operative use, or sterile for use in the operating room
- Surgical prep resistant
- Made of lasting gentian violet ink that is non-toxic and non-tattooing
- Latex free



PRODUCT SELECTION


ITEM #	TYPE	DESCRIPTION	PACKAGING
S-3000	Sterile	Multi-Tip fine marker with ruler and label set	50 / Box
S-3001	Sterile	Multi-Tip fine marker	50 / Box
S-30005	Sterile	Multi-Tip fine marker with ruler and label set (Volume Packaging)	500 / Box
S-30015	Sterile	Multi-Tip fine marker (Volume Packaging)	500 / Box
S-93000	Non-Sterile	Multi-Tip fine marker with ruler and label set	1000 / Box
S-93001	Non-Sterile	Multi-Tip fine marker	1000 / Box
S-94000	Non-Sterile	Ruler	100 / Bag
S-94001	Non-Sterile	Ruler (Bulk)	20,000 / Box
S-95000	Non-Sterile	Set of 3 labels	100 / Bag

SAVE by ordering our SKINS™ MARKERS in volume packaging. Cost effective and reduces waste.

BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



PRODUCT SELECTION

ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS + SLEEVE	ITEM # PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS
14-1028	914-1028 - SLV	914-1028	Yellow	0.028" / 0.71 mm
14-1035	914-1035 - SLV	914-1035	Light Blue	0.035" / 0.89 mm
14-1045	914-1045 - SLV	914-1045	White	0.045" / 1.14 mm
14-1054	914-1054 - SLV	914-1054	Dark Red	0.054" / 1.40 mm
14-1062	914-1062 - SLV	914-1062	Green	0.062" / 1.57 mm
14-1078	914-1078 - SLV	914-1078	Pink	0.078" / 1.98 mm
14-1093	914-1093 - SLV	914-1093	Grey	0.093" / 2.36 mm
14-1109	914-1109 - SLV	914-1109	Beige	0.109" / 2.77 mm
14-1125	914-1125 - SLV	914-1125	Purple	0.125" / 3.17 mm
14-1141	914-1141- SLV	914-1141	Brown	0.141" / 3.58 mm
14-1156	914-1156 - SLV	914-1156	Dark Blue	0.156" / 3.96 mm
14-1188	914-1188 - SLV	914-1188	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box

BATRIK SAFE-SOAK™ SPAC

Our Safe-Soak™ SPAC floor suction disk quickly removes fluid waste from the operating room floor through the application of suction.



- Highly durable: it won't be damaged if stepped on
- Easily maneuverable by foot around operating room floor
- Latex free, disposable, and non-sterile
- The bottom opening permits fluid collection through a highly absorbent material that comes into contact with its surface
- The hose is securely affixed on its side to reduce risk of catching equipment
- Its highly absorbent base material allows for quick fluid wicking and removal when suction is applied

BATRIK SAFE-SOAK™ ABSORBENT FLOOR MATS

Our Safe-Soak™ Absorbent Floor Mats are designed to absorb a variety of fluids from the operating room floor quickly and efficiently, ensuring the floor beneath remains slip-free, safe, and dry during surgical procedures.

- Highly absorbent
- Facilitates quick and easy clean up
- Latex free, single-use, disposable, and non-sterile
- Conforms to uneven surfaces
- Skid-resistant backing (blue absorbent mat only)




PRODUCT SELECTION

ITEM #	DESCRIPTION	ABSORBENCY	SIZE	PACKAGING
SS-20-1000	Safe-Soak™ / White extra absorbent mat	4.5 L / mat	28" x 48" (71 cm x 122 cm)	20 / Box
SS-20-2000	Safe-Soak™ / Blue absorbent mat with skid-resistant backer	2.5 L / mat	28" x 40" (71 cm x 102 cm)	20 / Box
SS-20-3000	Safe-Soak™ / Yellow absorbent mat	3 L / mat	30" x 52" (76 cm x 132 cm)	15 / Box
SS-21-1000	Safe-Soak™ SPAC	N/A	8.625" (22 cm diameter)	12 / Box

BATRIK MAGNETIC MAT

Our Magnetic Mat provides a sterile, neutral field for the safe placement and hands-free transfer of instruments. It is latex free, reusable, and autoclavable (up to 134° C).



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
MAG-02-1	Magnetic Mat	12" x 15" / 30 cm x 38 cm	1 Unit
MAG-02-5	Magnetic Mat (Volume Packaging)	12" x 15" / 30 cm x 38 cm	5 Units

BATRIK TOURNY™ + TOURNIQUETS

Our Tourniquets are comfortable and easy to use on all patients, including the elderly, undernourished, children, and newborns.



AVAILABLE COLORS



- Available in different sizes including compact sizes, which are convenient for phlebotomy trays
- Innovative latex free material has a more comfortable texture for patients, puts less stress on skin, and safeguards hospital staff and patients from latex allergens

SINGLE-USE BENEFITS

- Reduced risk of cross-contamination (patient to patient or patient to healthcare provider)
- Reduced risk of spread of nosocomial infections
- May improve your facility's compliance with OSHA's 1910-1030 (d)(1) "Blood-borne Pathogens Standard"

FOR YOUR CONVENIENCE:
25 UNITS PER ROLL
COME IN A DISPENSER BOX
FOR EASY ACCESS



PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	DIMENSIONS	PACKAGING
K-002	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 / Bag
K-004	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	100 / Box
K-009	TOURNY™+ Non-Sterile	Blue	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-010	TOURNY™+ Non-Sterile	White	0.200" x 22" / 0.5 cm x 56 cm	20 / Box
K-400	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-400PB	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	16 rolls 400 / Box
K-100	TOURNY™+ Non-Sterile	Assorted (3 rolls of each) Blue, Orange, Green, and Pink	1" x 18" / 2.5 cm x 45 cm	12 rolls 300 / Box
K-254	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	250 / Box
K-2-18RL	TOURNY™+ Non-Sterile	Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-18PBRL	TOURNY™+ Non-Sterile (Pebbled)	Light Blue	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-3-18RL	TOURNY™+ Non-Sterile	Green	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-6-18RL	TOURNY™+ Non-Sterile	Pink	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-7-18RL	TOURNY™+ Non-Sterile	Orange	1" x 18" / 2.5 cm x 45 cm	1 Roll
K-2-29RL	Long TOURNY™+ Non-Sterile	Blue	1" x 29.5" / 2.5 cm x 75 cm	1 Roll

BATRIK
SAVEFACE™ SHOE COVERS

- Our SaveFace™ superior medical shoe covers are durable, safe, and easy to use. Made of polypropylene fabric with robust welds, our shoe covers offer the ultimate in durability. The triple-elastic run and anti-skid print allows for a smooth fit around the foot and reduces the risk of slipping.
- Particle count tested and complies with IEST-RP-CC003.4
- 100% latex and silicone free


 MADE IN CANADA 
PRODUCT SELECTION

ITEM #	SIZE	PACKAGING
SF-SC02UV	Universal Up to shoe size 11	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
SF-SC02UV-8	Universal Up to shoe size 11	Volume Packaging 100 / bag - 8 bags / case
SF-SC02XL	Extra-Large Size 11 - 15	Regular Packaging 200 / dispenser box - 4 dispenser boxes / case
SF-SC02XL-8	Extra-Large Size 11 - 15	Volume Packaging 100 / bag - 8 bags / case

BATRIK WTFOG™ by G.O.L.F.F.™ REUSABLE ANTI-FOG CLOTHS

“WHAT THE FOG?!”

We've all experienced it.

Eyewear fogging up while wearing a mask is not only unpleasant and highly frustrating, but it can be dangerous, too.

INTRODUCING THE FOG FIX:

Our WTFog™ reusable anti-fog cloths prevent fogging on all types of eyewear. Ideal for use when wearing a mask or face shield, WTFog™ takes you from “What The Fog?!” to confidently ready for your next patient in a snap.

WHY YOU'LL LOVE IT:

- Eliminates fogging on all glass, polycarbonate, and plastic eyewear lenses
- Each reusable nano-infused cloth lasts up to 600 wipes or 6 months
- Microfiber cloths are infused with our proprietary heavy duty anti-fog treatment
- Quick and easy application - just wipe and wear
- Lint-free, scratch-free, no greasy smearing
- Long lasting protective coating provides effective protection against fog, humidity, sweat, and fast temperature changes.
- 100% Natural Titanium Dioxide (TiO₂) Nano Formula developed in Germany meets strict REACH standards
- Non-Sterile



“To be able to complete a whole appointment without adjusting my eye protection makes maintaining COVID PPE protocols that much easier.”

DR. MOLLY FRIEDMAN

HOLLAND BLOORVIEW KIDS REHABILITATION HOSPITAL
 TORONTO, ON



PERFECT WHEN:

- You wear glasses/goggles and a mask/face shield
- Your profession regularly exposes you to fogging and you need crystal clear vision without the risk or hassle of fogging up



DIRECTIONS FOR USE:

Open resealable pouch and remove the WTFog™ reusable anti-fog cloth. Wipe cloth on clean surface (glasses/goggle lens) being careful not to remove the protective coating delivered by the cloth. There is no wait time required before wearing your protected item. Enjoy your fog-free experience!

CLOTH STORAGE & CARE:

After use, keep cloth in the resealable pouch until next use, ensuring the seal is fully closed. If cloth becomes too dry, add up to 1 ml of water to pouch to reactivate cloth. For best results, store at ambient temperature/humidity. Dispose cloth after 600 uses or 6 months.

**External surface use only.
Keep out of reach of children.
Do not wash**

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
WTFOG-30	30 WTFog™ in a dispenser box measuring: 3 7/16" x 5 2/16" x 6 2/16" (8.7 cm x 13 cm x 15.4 cm)	30 WTFog™ / dispenser box
WTFOG-100	100 WTFog™ - bulk packed	100 WTFog™ / box
WTFOG-1200	1200 WTFog™ - bulk packed	1200 WTFog™ / case (12 boxes of 100 bulk packed)

1 WTFog™ Cloth per Resealable Pouch
 WTFog™ Cloth: 5.5" x 5.5" (13.9 cm x 13.9 cm)
 WTFog™ Resealable Pouch: 3.5" x 5" (8.9 cm x 12.7 cm)

IMPERIAL SURGICAL LTD. STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, Imperial Surgical Ltd.'s patented Stackable Step Stools have been designed with infection prevention as a prime focus.



PATENTED

- Weight of stool: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Antimicrobial coating prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast, and mildew
- Cart washable. Model OR-36364 (without mat) can be sterilized
- Anti-slip surface for enhanced stability
- Rubber-free corners avoids fluid collection
- No welds or large hollow openings to trap debris or fluids
- Transport cart available
- 300 lb (136 kg) weight capacity
- 7 year warranty covers all parts (step stool only)

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR
OR-36363	Stackable Step Stool with mat	Black
OR-36363-02	Stackable Step Stool with mat	Dark Blue
OR-36363-03	Stackable Step Stool with mat	Green
OR-36363-07	Stackable Step Stool with mat	Orange
OR-36364	Stackable Step Stool without mat	N/A

PACKAGING: Sold by Unit (Dealers: 4 / Box)

OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING
OR-3693	Transport cart	1 Unit
809-1705	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box
899-1710	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box

REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING
OR-36363-1A	Replacement mat for OR-36363	1 Unit
OR-36363-1B	Replacement Foot Pegs	1 Unit

TIP: COLOR-CODE YOUR STEP STOOLS BY DEPARTMENT



BATRIK BLUE CLOUD-GEL™ ANTI-FATIGUE MATS


Our Blue Cloud-Gel™ Anti-Fatigue Mats have been scientifically engineered for use in the O.R. using patented gel technology and premium foam to provide the ultimate comfort and support for surgeons, nurses, and technicians particularly when movement is restricted for prolonged durations.

- Proprietary 3/4" Blue Cloud-Gel™ foam technology will not "bottom-out" like conventional foam mats
- Conforms to and supports feet promoting efficient circulation, reducing physical discomfort and fatigue
- Promotes proper posture while diminishing back, knee, and foot discomfort
- Latex free, water resistant, and nonabsorbent
- Reusable, easy to clean, and Cart Washer compatible (contact us for more details)
- Excellent friction properties when damp or dry surgical booties are worn
- Completely sealed with an RF weld for excellent protection against moisture absorption
- US Standards: CALL 117 and NFPA 260 Compliant
- 1 Year Warranty

DID YOU KNOW? STUDIES HAVE DEMONSTRATED THAT ANTI-FATIGUE MATS IMPROVE PERFORMANCE, ENERGY, AND POST-OPERATIVE DISCOMFORT.

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
800-2401	18" X 24" / 46 cm x 61 cm	1 / Box
800-2405	18" X 24" / 46 cm x 61 cm	5 / Box
800-3201	20" X 32" / 51 cm x 81 cm	1 / Box
800-3203	20" X 32" / 51 cm x 81 cm	3 / Box
800-4801	20" X 48" / 51 cm x 122 cm	1 / Box
800-4802	20" X 48" / 51 cm x 122 cm	2 / Box
800-7201	30" X 72" / 76 cm x 183 cm (Special Order)	1 / Box




BATRIK BLUE-CLOUD NIMBUS™ ANTI-FATIGUE MATS

Our Blue-Cloud Nimbus™ Anti-Fatigue Mats are designed for general hospital use. Made with proprietary high-density “energy-return” polyurethane materials, they provide unmatched comfort, relief, and ergonomic support.

- Our Energy-Return properties will help reduce leg, foot, lower back discomfort, and work-related injuries
- Longer lasting support as 3/4” proprietary foam blend will not bottom out like most other foam based anti-fatigue mats
- Antimicrobial properties blended into the mat provides added protection
- Highly water resistant surface will not absorb water or dirt
- Lightweight and easy to clean
- Textured design on surface and high-traction base ensures greater safety. (US Standard NFSI 101-C)
- Beveled edging reduces tripping hazards and anti-slip backing enhances safety
- Durable design is ideal for commercial applications, yet still ergonomically and aesthetically pleasing for home or office environments
- Flammability/Burn resistant per CFR 16 Parts 1630/1631
- Made with eco-friendly materials, is non-toxic, and latex free
- 5 Year Warranty (non-industrial usage), 4 Year Warranty (industrial usage)

GOOD TO KNOW: IF WITHIN 30 DAYS YOU ARE NOT 100% SATISFIED, CONTACT YOUR DEALER FOR A REFUND (EXCLUDES TRANSPORT, HANDLING, AND TAXES)



Available in Taupe or Black

PRODUCT SELECTION

ITEM #	DIMENSIONS INCHES / CM	PACKAGING
BLACK COLOR		
809-3201	20" X 32" / 51 cm x 81 cm	1 / Box
809-3204	20" X 32" / 51 cm x 81 cm	4 / Box
809-4801	20" X 48" / 51 cm x 122 cm	1 / Box
809-4804	20" X 48" / 51 cm x 122 cm	4 / Box
809-3601	24" X 36" / 61 cm x 91 cm	1 / Box
809-7201	20" X 72" / 51 cm x 183 cm	1 / Box
809-6001	36" X 60" / 91 cm x 152 cm	1 / Box
TAUPE COLOR		
810-3201	20" X 32" / 51 cm x 81 cm	1 / Box
810-3204	20" X 32" / 51 cm x 81 cm	4 / Box
810-4801	20" X 48" / 51 cm x 122 cm	1 / Box
810-4804	20" X 48" / 51 cm x 122 cm	4 / Box
810-3601	24" X 36" / 61 cm x 91 cm	1 / Box
810-7201	20" X 72" / 51 cm x 183 cm	1 / Box
810-6001	36" X 60" / 91 cm x 152 cm	1 / Box

BATRIK G.O.L.F.F.™ ANTI-FOG SOLUTION

Our G.O.L.F.F.™ Anti-Fog Solution solves the fogging problem often encountered during Laparoscopic surgical procedures.

- Prevents the occurrence of fog on lenses
- Ensures clearer vision and avoids unnecessary interruptions during surgery
- Set includes an easy-squeeze 6 mL bottle and an adhesive radiopaque sponge for secure placement
- OEM (available)
- Sterile
- Latex free



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
G.O.L.F.F.	Anti-Fog Solution & Adhesive Sponge set	20 / Box


BATRIK ENDOSCOPY BO2™ BITE BLOCKS WITH NASAL OXYGEN DELIVERY

Our BO2™ Bite Blocks significantly enhance patient safety by providing supplemental oxygenation during gastrointestinal endoscopy, preventing hypoxia in sedated patients, and ensuring superior patient comfort.

- Saves time and money by eliminating the need for separate nasal prongs and oxygen masks
- Designed to be free of sharp edges to ensure patient comfort and safety
- Available exclusively in Single-Use
- Optional retention strap is available to ensure a secure fit
- The Bite Blocks and optional retention straps are all latex free, lead free and DEHP free

BO2
(FRONT SIDE)

BO2
(BACK SIDE)



PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
BO2	Single-Use	18mm / 54FR, includes strap	100 / Box
BO2-N	Single-Use	18mm / 54FR, without strap	100 / Box

BATRIK ENDOSCOPY BITESTRONG™ BITE BLOCKS

ADULT & PEDIATRIC BITE BLOCKS

Our Bitestrong™ Bite Blocks provide a safe passage through the patient's mouth during TEE, protecting the patient's teeth and endoscope.

Our molding process ensures a smooth finish with no sharp edges (flashing) that can often irritate a patient's mouth. All Bite Blocks are individually wrapped in sealed, easy to open PE bags and are latex free, lead free, and DEHP free. Additional straps are available, if required

SINGLE USE

- Available without dental retention

REUSABLE

- Designed for repeated usage and will withstand commonly practiced chemical and steam sterilization cycles (134°C / 275°F).
- Available with or without an adjustable retention strap

PRODUCT SELECTION

ITEM #	TYPE	DESCRIPTION	PACKAGING
B1	Single-Use Pediatric	16 mm / 48FR, includes strap	100 / Box
B1-25	Single-Use Pediatric	16 mm / 48FR, includes strap	25 / Box
B1-N	Single-Use Pediatric	16 mm / 48FR, without strap	100 / Box
B3	Single-Use	20 mm / 60FR, includes strap	100 / Box
B3-N	Single-Use	20 mm / 60FR, without strap	100 / Box
B3-SA	Single-Use	20 mm / 60FR, with strap attached to one side	100 / Box
B3-NDR	Single-Use	20 mm / 60FR, no dental retention rim	100 / Box
B3-NDR-N	Single-Use	Same as B3-N, no dental retention rim, without strap.	100 / Box
B4	Single-Use	19 mm / 60FR	100 / Box
B4T	Single-Use	19 mm / 60FR, includes strap	100 / Box
B4T-N	Single-Use	19 mm / 60FR, without strap	100 / Box
B6	Reusable	19 mm / 60FR	100 / Box
B6T	Reusable	19 mm / 60FR, includes strap	100 / Box
B6T-N	Reusable	19 mm / 60FR, without strap	100 / Box
B6-25	Reusable	19 mm / 60FR	25 / Box
B6T-25	Reusable	19 mm / 60FR, includes strap	25 / Box
B6T-N-25	Reusable	19 mm / 60FR, without strap	25 / Box
BS	Single-Use	Bite-Block straps 21.7" x .79" / 55 cm x 2 cm	100 / Box



B5



B1



B3



B4



B6

BATRIK COUVRETTE™ VALVE COVERS

Our Couvrette™ Valve Covers eliminate the need for manual cleaning and reprocessing of reusable valves. Our single use Valve Covers effectively reduce cross contamination between patients and also reduce fluid leakage from the parts during procedures.

- Easily and securely fit to all valves
- Reduce the potential for errors
- Single-Use
- Non-Sterile
- Latex free



16-0044

TIP: REUSABLE VALVE COVERS ARE NOTORIOUSLY DIFFICULT TO CLEAN. SAVE TIME AND REDUCE RISK WITH COUVRETTE™.



16-0022

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR	PACKAGING
16-0022	Biopsy Valve Cover - Olympus® & Fujinon®	Blue	100 / Box
16-0044	Biopsy Valve Cover - Pentax®	Pink	100 / Box


NOTE: Contact us to verify compatibility with your scope

BATRIK DAYBRUSH™ CHANNEL BRUSHES

REUSABLE

Our Reusable DayBrush™ Channel Brushes can be used for cleaning endoscopes for an entire day. In fact, it's been proven to effectively clean up to 10 scopes before the need to discard.*

- Soft tip ensures easy manipulation within the channel with no risk of damage
- The back-end is sealed and softened to prevent fluid from getting trapped inside the tube. The softening process also prevents damage if the cleaning brush is used in reverse.
- Can be steam sterilized.
- Bristle volume is 20mm in length – nearly double most brushes on the market. A large brush provides more surface contact during the cleaning process.
- Significant per use savings when compared to other reusable or disposable channel cleaning brushes.



() Discard after cleaning 10 scopes or sooner if any brush degradation is noticed. After ten applications, the brush may have the appearance that it can continue to clean effectively, but we discourage prolonged use due to risk of cross contamination. Actual number of uses will vary based on handling and processing.*

PRODUCT SELECTION

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
30-18050-T	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-T	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Single
30-23050-2T	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.07" / 1.75 mm	Double

PACKAGING: 50 / Box

OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.


"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

BATRIK ARC™ CHANNEL BRUSHES SINGLE-USE

BATRIK offers one of the largest selection of Endoscopy Channel Cleaning Brushes in the industry. Our brushes are available in a variety of lengths, materials, and styles to accommodate the widest range of endoscopic devices used in the healthcare industry today.

No matter what your requirements, you'll be able to find the brush that suits your specific needs. If you don't see the brush you require, simply call us. We have over 50 different brushes available and can custom fabricate to your precise specifications.



OPTIONAL KIT ADD-ONS

ADD THE SUFFIX:

"F" to receive a kit with our Flat Sponge (Packaging: Box of 40). Applicable on certain models only, call for details.

"K" to receive a kit with the single use Valve / Control Cleaning Brush 5/11 Double Brushes (item # 27-57787). Applicable on all models.

"FK" to receive both the sponge and valve cleaning brushes.

SERIES 333 & 311



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
311-23050	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-23050-2	90" / 230 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Double
311-18050	70" / 180 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single
311-10050	39" / 100 cm	0.94" / 24 mm	0.20" / 5 mm	0.073" / 1.85 mm	Single

PACKAGING: 50 / BOX

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.85 MM - PUNCHY PINK POLYPROPYLENE

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
333-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Single
333-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.073" / 1.85 mm	Double

PACKAGING: 50 / Box

SERIES 300 & 322



PRODUCT SELECTION

5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
300-18050	70" / 180 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
300-23050	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Single
300-23050-2	90" / 230 cm	0.79" / 20 mm	0.20" / 5 mm	0.070" / 1.75 mm	Double


PACKAGING: 50 / Box

5/6 MM TAPERED BRUSHES – SHEATH Ø 1.75 MM - METALLIC GREY DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
322-23056	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Single
322-23056-2	90" / 230 cm	0.79" / 20 mm	0.20", 0.24" / 5, 6 mm	0.070" / 1.75 mm	Double

PACKAGING: 50 / Box

SERIES 344




PRODUCT SELECTION

3 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.55 MM - METALLIC PURPLE DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
344-23030	90" / 230 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
344-18030	70" / 180 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single
344-10030	39" / 100 cm	0.79" / 20 mm	0.12" / 3 mm	0.060" / 1.55 mm	Single

PACKAGING: 50 / Box

SERIES 35



PRODUCT SELECTION

2.5 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-23025	90" / 230 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single
35-12025	47.2" / 120 cm	0.59" / 15 mm	0.098" / 2.5 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

2 MM CYLINDRICAL BRUSHES – SHEATH Ø 1.00 MM - YELLOW DELRIN

ITEM #	LENGTH	BRUSH LENGTH	BRUSH (Ø)	SHEATH (Ø)	BRUSH ENDS
35-23020	90" / 230 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-12020	47.2" / 120 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single
35-18020	70.8" / 180 cm	0.59" / 15 mm	0.079" / 2.0 mm	0.039" / 1.00 mm	Single

PACKAGING: 50 / Box

BATRIK VALVE / CONTROL BRUSHES


27-05197



27-57787



40-18511



36-10010



27-15311




PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH	BRUSH 1 LENGTH	BRUSH 2 LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)
27-05197	Double Trumpet 5/7 White Color	5.00" / 13 cm	0.20" / 5 mm	1.00" / 25 mm	0.20" / 5 mm	0.28" / 7 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-05197-50						
27-57787	Double. Universal 5/11 Yellow Color	6.00" / 15 cm	1.18" / 30 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-57787-50						
40-18511	Double Universal 5/11 Extra Long Blue Color	8.07" / 20.5 cm	1.57" / 40 mm	0.79" / 20 mm	0.43" / 11 mm	0.20" / 5 mm
PACKAGING: Bag of 3 (reusable and can be steam sterilized). Volume Packaging: Bag of 50 use ITEM # 40-18511-50						
36-10010	Single end large Red Color	4.3" / 11 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 36-10010-50						
27-15311	Single end large Blue Color	4" / 10.3 cm	1.57" / 40 mm	N/A	0.39" / 10 mm	N/A
PACKAGING: Bag of 3 (single-use recommended). Volume Packaging: Bag of 50 use ITEM # 27-15311-50						

BATRIK V-BRUSHES EXTERNAL CLEANING BRUSHES

Our External Cleaning Brushes for endoscopes / optical instruments provide excellent cleaning capabilities. With their unique U-shape bristles and innovative design, our brushes are the ultimate tool for cleaning the exterior of your endoscopic instruments. Free from sharp or metallic ends that can scratch or cause damage to valuable equipment, both single-use and reusable models are latex free and feature ergonomic handles.





PRODUCT SELECTION

ITEM #	DESCRIPTION	BRUSH LENGTH / (Ø)	HANDLE LENGTH / (Ø)	TOTAL LENGTH	PACKAGING
36-20005	Single-Use External Cleaning Brush	70 mm / 25 mm	80 mm / 3 mm	150 mm	50 / Box
36-20010	Reusable External Cleaning Brush	90 mm / 25 mm	120 mm / 12 mm	210 mm	3 / Bag


BATRIK SILICONE PIN MAT

Our Silicone Pin Mat cushions, protects, and organizes delicate surgical instruments for sterilization. The elevated nubs provide greater air-flow ventilation and drying, significantly reducing residual moisture collection.

- Reusable
- Autoclavable
- Latex free



FRONT SIDE
Pins keep instruments securely in place



BACK SIDE
Numerous holes and elevated nubs provide superior venting and reduce water collection

PRODUCT SELECTION

ITEM #	DIMENSIONS	PACKAGING
13-2001	13" x 18" (30 cm x 45 cm)	1 unit
13-2005	13" x 18" (33 cm x 45 cm)	5 units

GOOD TO KNOW:
CUSTOM SIZES
ARE AVAILABLE
UPON REQUEST

BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




NEW!

18-CP06

COLOUR	BLUE	RED	PINK
USAGE	SINGLE-USE	REUSABLE <i>withstands repeated sterilization up to 134°C</i>	REUSABLE <i>withstands repeated sterilization up to 134°C</i>
DIMENSIONS	2.8" x 4" (7.3 cm x 10.3 cm)	2.8" x 4" (7.3 cm x 10.3 cm)	3.4" x 4.7" (8.5 cm x 12 cm)
LATEX-FREE	✓	✓	✓
NO MOISTURE RETENTION	✓	✓	✓
FABRICATION	Flexible, silicone-like plastic	Flexible, silicone-like plastic	Our most flexible plastic for easiest placement and most secure protection

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors feature foam that is 0.25" thicker than competitor products. They are especially useful with heavier instrument sets.



18-CPF-02

PRODUCT SELECTION



ITEM #	DESCRIPTION	PACKAGING
18-CP02**	Corner Pockets™ Single-Use, Blue	240 / Bag
18-EVAL2**	Corner Pockets™ Single-Use, Blue	24 / Bag
18-CP07	Corner Pockets™ Reusable, Red	80 / Bag
18-CORNER	Corner Pockets™ Reusable, Red	320 / Box (4 Bags of 80 pcs)
18-EVAL7	Corner Pockets™ Reusable, Red	24 / Bag
18-CP06	Corner Pockets™ Reusable, Pink	80 / Bag
18-CPF-02	Soft Pockets™ Single-Use Foam Corner Protectors	80 / Box (24 Boxes / Pallet)

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.

BATRIK INSTRUMENT PROTECTORS

Our Instrument Protectors protect the tips and edges of instruments to prevent damage during handling and sterilization. Available in a wide array of color coded sizes for almost every application, including osteotome, retractor, double, and perforated.

- Steam and gas permeable
- Perforated protectors reduce vacuum seal and moisture build-up
- Validated for pre-vacuum sterilization
- Multiple color-coded sizes
- Latex Free
- Non-Sterile



PRODUCT SELECTION

REGULAR & PERFORATED / CLEAR



RC-07



PC-07

ITEM # REGULAR	ITEM # PERFORATED	COLOR	DIMENSIONS
RC-01	PC-01	Clear	0.062" x 0.75" (1.57 mm x 19 mm)
RC-02	PC-02	Clear	0.075" x 0.75" (2 mm x 19 mm)
RC-03	PC-03	Clear	0.109" x 0.75" (2.77 mm x 19 mm)
RC-04	PC-04	Clear	0.125" x 1" (3.18 mm x 25.4 mm)
RC-05	PC-05	Clear	0.187" x 1" (4.75 mm x 25.4 mm)
RC-06	PC-06	Clear	0.406" x 0.75" (10.3 mm x 19 mm)
RC-07	PC-07	Clear	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RC-08	PC-08	Clear	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RC-09	PC-09	Clear	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RC-10	PC-10	Clear	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / TINTED



RT-07



PT-07

RT-01	PT-01	Tinted White	0.062" x 0.75" (1.57 mm x 19 mm)
RT-02	PT-02	Tinted Blue	0.075" x 0.75" (2 mm x 19 mm)
RT-03	PT-03	Tinted Green	0.109" x 0.75" (2.77 mm x 19 mm)
RT-04	PT-04	Tinted Red	0.125" x 1" (3.18 mm x 25.4 mm)
RT-05	PT-05	Tinted Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RT-06	PT-06	Tinted Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RT-07	PT-07	Tinted Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RT-08	PT-08	Tinted Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RT-09	PT-09	Tinted Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RT-10	PT-10	Tinted Assorted	Assorted

PACKAGING: Bag of 100

REGULAR & PERFORATED / SOLID



RS-07



PS-07


RS-01	PS-01	Solid White	0.062" x 0.75" (1.57 mm x 19 mm)
RS-02	PS-02	Solid Blue	0.075" x 0.75" (2 mm x 19 mm)
RS-03	PS-03	Solid Green	0.109" x 0.75" (2.77 mm x 19 mm)
RS-04	PS-04	Solid Red	0.125" x 1" (3.18 mm x 25.4 mm)
RS-05	PS-05	Solid Yellow	0.187" x 1" (4.75 mm x 25.4 mm)
RS-06	PS-06	Solid Maroon	0.406" x 0.75" (10.3 mm x 19 mm)
RS-07	PS-07	Solid Orange	0.062" x 0.375" x 1" (1.57 mm x 9.53 mm x 25.4 mm)
RS-08	PS-08	Solid Brown	0.078" x 0.625" x 1" (2 mm x 15.875 mm x 25.4 mm)
RS-09	PS-09	Solid Black	0.125" x 1" x 1" (3.18 mm x 25.4 mm x 25.4 mm)
RS-10	PS-10	Solid Assorted	Assorted

PACKAGING: Bag of 100

ITEM #	DIMENSIONS
OP-11	0.25" x 0.75" (6.35 mm x 19.05 mm)
OP-12	0.375" x 1" (9.53 mm x 25.4 mm)
OP-13	0.5" x 1" (12.7 mm x 25.4 mm)
OP-14	0.75" x 1" (19.05 mm x 25.4 mm)
OP-15	1" x 1" (25.4 mm x 25.4 mm)
OP-16	1.25" x 1.25" (31.75 mm x 31.75 mm)
OP-17	1.5" x 1.5" (38.1 mm x 38.1 mm)
OP-18	1.75" x 1.25" (44.45 mm x 31.75 mm)
OP-20	Assorted


PACKAGING: Bag of 50

OSTEOTOME PROTECTORS



ITEM #	DIMENSIONS	PACKAGING
RGS-01	1.5" (38 mm)	55 / Bag
RGS-02	2" (51 mm)	45 / Bag
RGS-03	3" (76.2 mm)	40 / Bag


SINGLE-TIP RETRACTOR PROTECTORS



ITEM #	COLOR	SIZE
RGM-01	Pink	Small
RGM-02	Clear	Small
RGM-03	Green	Medium
RGM-04	Clear	Medium

PACKAGING: Bag of 50


WEITLANER RETRACTOR PROTECTORS



ITEM #	COLOR	DIMENSIONS
DG-01	Clear	0.0625" x 0.75" (1.59 mm x 19.05 mm)
DG-02	Blue	0.078" x 0.75" (1.98 mm x 19.05 mm)
DG-03	Green	0.109" x 0.75" (2.78 mm x 19.05 mm)
DG-04	Red	0.125" x 1" (3.18 mm x 25.4 mm)
DG-05	Yellow	0.187" x 1" (4.76 mm x 25.4 mm)
DG-07	Orange	0.062" x 0.375" x 1" (1.59 mm x 9.53 mm x 25.4 mm)
DG-10	Assorted	Assorted

PACKAGING: Bag of 50

PERFORATED DOUBLE PROTECTORS



BATRIK
PrAmLi Power³ LUBE™
Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. **PrAmLi Power³ LUBE™** is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:
 FOR MILK BATH SOLUTION,
 DILUTE 1 PART SOLUTION
 IN 6 PARTS WATER

Formulated with
PrAmLi Power³ Enzymes™




PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles

BATRIK S.T.I.C.S.

Our Simple Tracking Identification Color System, or S.T.I.C.S., is an easy to use, latex free, tool that allows you to quickly identify your instruments by department and/or procedure. The tape adheres securely and every tape sheet comes with a protective cover to ensure it lasts as long as you need.



Each 8.5" x 11" (216 mm x 279 mm) sheet contains the following adhesive strips:

- 1" (25.4 mm) x 220
- 1.5" (38.1 mm) x 44
- 2" (50.8 mm) x 44
- Total tape per sheet 374" (9500 mm)

PRODUCT SELECTION

SOLID COLORS	
ITEM #	COLOR
T-850-2-01	White
T-850-2-02	Royal Blue
T-850-2-03	Kelly Green
T-850-2-04	Red
T-850-2-05	Yellow
T-850-2-06	Purple
T-850-2-07	Orange
T-850-2-08	Brown
T-850-2-09	Black
T-850-2-10	Gray
T-850-2-11	Lavender
T-850-2-12	Beige
T-850-2-13	Light Blue
T-850-2-14	Lime Green
T-850-2-15	Pink
T-850-2-17	Plum
T-850-2-18	Mint Green

BLACK STRIPES	
ITEM #	COLOR
T-850-3-01	Black/White
T-850-3-03	Black/Kelly Green
T-850-3-04	Black/Red
T-850-3-05	Black/Yellow
T-850-3-07	Black/Orange
T-850-3-10	Black/Gray
T-850-3-11	Black/Lavender
T-850-3-12	Black/Beige
T-850-3-13	Black/Light Blue
T-850-3-14	Black/Lime Green
T-850-3-15	Black/Pink
T-850-3-17	Black/Plum

TRI-COLOR	
ITEM #	COLOR
T-850-5-51	Red/White/Blue
T-850-5-54	Red/Black/Green

WHITE STRIPES	
ITEM #	COLOR
T-850-4-02	White/Royal Blue
T-850-4-03	White/Kelly Green
T-850-4-04	White/Red
T-850-4-05	White/Yellow
T-850-4-07	White/Orange
T-850-4-08	White/Brown
T-850-4-11	White/Lavender
T-850-4-13	White/Light Blue
T-850-4-14	White/Lime Green
T-850-4-15	White/Pink
T-850-4-17	White/Plum

FOUR-COLOR	
ITEM #	COLOR
T-850-2-31	Red/Yellow/ Purple/Orange
T-850-2-32	Blue/White/ Black/Green


SPECIALTY	
ITEM #	COLOR
T-850-7-04	Yellow/Red

BATRIK ID TAPE



A simple solution for marking surgical instruments with the aid of color-coding, our ID Tape is durable and will not chip, flake, or peel. It can be used in autoclaves, withstanding repeated exposure to sterilization. Independent validation tests demonstrated the efficiency of sterilization of marked instruments, even underneath the Batrik ID Tape.

AVAILABLE WIDTHS



Quickly identifies instruments by department or procedure



Easy to apply



Withstands repeated exposure to sterilization



Has more per tape roll (300" / 7620 mm)



Available in 9 colors and 2 widths



Conveniently packaged in a clear box



Non-Sterile

PRODUCT SELECTION

ITEM #	ITEM#	COLOR
0.125" / 318 MM	0.25" / 635 MM	
T-125-01	T-250-01	White
T-125-02	T-250-02	Blue
T-125-03	T-250-03	Green
T-125-04	T-250-04	Red
T-125-05	T-250-05	Yellow
T-125-06	T-250-06	Purple
T-125-07	T-250-07	Orange
T-125-08	T-250-08	Brown
T-125-09	T-250-09	Black

BATRIK TAPE-B-GONE™ TAPE REMOVER

Tape-B-Gone™ is a uniquely formulated and ready-to-use proprietary formula that easily removes adhesives from medical equipment and instruments. Unlike many other brands, Tape-B-Gone™ does not contain highly dangerous or caustic chemicals, is not flammable or combustible, and is safe for use on metal, glass, plastic, and rubber without any special equipment or ventilation.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
T-08RMV	8 oz (237 ml) bottle	1 Unit
T-32RMV	8 oz (237 ml) bottle	4 / Pack
T-SAMP-RMV	0.5 oz (15 ml) bottle	1 Unit

BATRIK PADLOCKS SECURITY TAGS

Our Padlock Security Tags for instrument containers and general hospital applications are specifically designed for the healthcare market. Our material selection, fabrication process, and engineering ensure consistent quality and high temperature compliance with the various steam sterilization temperatures in worldwide hospital applications.

Available with optional reactive ink indicator dots that provide visual sterility indication. Light blue dots turn brown/dark brown when steam sterilized.



Made of polypropylene



Latex and DEHP Free



Will not melt or soften with any established steam sterilization protocol




Fully open design (45°) facilitates insertion into instrument containers



Easy twist breaking design



Large base for easier handling




Available in different size




Custom Colors and OEM available upon request

AVAILABLE IN THREE FORMATS


700-SERIES



800-SERIES



900-SERIES



PRINTED & CUSTOM NUMBERED PADLOCKS

Available in 700, 800 & 900 SERIES

PRODUCT SELECTION

COLOR	PRINTED WORD
PRINTED PADLOCKS	
WHITE	CLEAN
BLUE	STERILE
RED	REPAIR
YELLOW	SHARPEN
YELLOW	LOANER
BLACK	DISCARD
CUSTOM NUMBERED	
ALL COLORS	up to 8 numbers



CUSTOM NUMBERED
Custom printed padlocks (up to 8 characters)



PRINTED
6 different printed labels available

NOTE: Printed & Custom Numbered Padlocks are available in 700, 800, and 900 Series and can be ordered in 100, 200, and 500 unit quantities. Please contact us for ordering details.

INDICATOR DOTS



REGULAR



700-SERIES PADLOCKS

The 700 series padlocks are perfect for "Aesculap" containers and those with a short, large opening.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701	15-711	15-710
BLUE	15-702	15-712	15-720
GREEN	15-703	15-713	15-730
RED	15-704	15-714	15-740
YELLOW	15-705	15-715	15-750
PINK	15-706	15-716	15-760
ORANGE	15-707	15-717	15-770
BLACK	15-709	15-719	15-790
LIGHT BLUE	15-722	15-732	15-742

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-701-IND	15-711-IND	15-710-IND
BLUE	15-702-IND	15-712-IND	15-720-IND
GREEN	15-703-IND	15-713-IND	15-730-IND
RED	15-704-IND	15-714-IND	15-740-IND
YELLOW	15-705-IND	15-715-IND	15-750-IND
PINK	15-706-IND	15-716-IND	15-760-IND
ORANGE	15-707-IND	15-717-IND	15-770-IND
BLACK	15-709-IND	15-719-IND	15-790-IND
LIGHT BLUE	15-722-IND	15-732-IND	15-742-IND

Double locking hooks can be inserted in a large guide hole



The grooves on the padlock facilitate handling



INDICATOR DOTS

REGULAR



800-SERIES PADLOCKS

The 800 series padlocks have a longer arm (32 mm) and have a single, thinner hook than the 700 series. They are perfect for almost any type of container, even those with long and thin channels.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801	15-811	15-810
BLUE	15-802	15-812	15-820
GREEN	15-803	15-813	15-830
RED	15-804	15-814	15-840
YELLOW	15-805	15-815	15-850
PINK	15-806	15-816	15-860
ORANGE	15-807	15-817	15-870
BLACK	15-809	15-819	15-890
LIGHT BLUE	15-822	15-832	15-842

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-801-IND	15-811-IND	15-810-IND
BLUE	15-802-IND	15-812-IND	15-820-IND
GREEN	15-803-IND	15-813-IND	15-830-IND
RED	15-804-IND	15-814-IND	15-840-IND
YELLOW	15-805-IND	15-815-IND	15-850-IND
PINK	15-806-IND	15-816-IND	15-860-IND
ORANGE	15-807-IND	15-817-IND	15-870-IND
BLACK	15-809-IND	15-819-IND	15-890-IND
LIGHT BLUE	15-822-IND	15-832-IND	15-842-IND

Single locking hook can be inserted in a large guide hole




The locking system allows for easy closure



INDICATOR DOTS

REGULAR



900-SERIES PADLOCKS

The 900 series padlocks are identical to the 800 series, but with a larger rectangular base which allows for the placement of a 20 mm x 30 mm tracking label.

PRODUCT SELECTION

REGULAR PADLOCKS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901	15-911	15-910
BLUE	15-902	15-912	15-920
GREEN	15-903	15-913	15-930
RED	15-904	15-914	15-940
YELLOW	15-905	15-915	15-950
PINK	15-906	15-916	15-960
ORANGE	15-907	15-917	15-970
BLACK	15-909	15-919	15-990
LIGHT BLUE	15-922	15-932	15-942

WITH INDICATORS DOTS

COLOR	100 UNITS	200 UNITS	500 UNITS
WHITE	15-901-IND	15-911-IND	15-910-IND
BLUE	15-902-IND	15-912-IND	15-920-IND
GREEN	15-903-IND	15-913-IND	15-930-IND
RED	15-904-IND	15-914-IND	15-940-IND
YELLOW	15-905-IND	15-915-IND	15-950-IND
PINK	15-906-IND	15-916-IND	15-960-IND
ORANGE	15-907-IND	15-917-IND	15-970-IND
BLACK	15-909-IND	15-919-IND	15-990-IND
LIGHT BLUE	15-922-IND	15-932-IND	15-942-IND

Single locking hook
for use with smaller diameters



Ideal for adding labels/ indicators indicators or markings to the large rectangular base





BATRIK SAVEFACE™ FACE SHIELD

Our SaveFace™ Face Shield (Visor) is lightweight and allows for good peripheral vision while acting as a physical barrier against splashes. Can be used for extended periods due to comfort-focused design and superior quality materials.

Choose amongst two models:

- The classic **non-draped** model closes the gap between the shield and the chest bone when head is tilted forward for added protection.
- **Draped** offers full-length coverage extending from the top of the eyebrow to below the chest, covering the anterior neck area as per CDC recommendations.



ASSEMBLED IN CANADA

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-100	SAVEFACE™ no drape	Shield Length: 9" / 23 cm	9" (23 cm)	100 / Box
SF-180	SAVEFACE™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm Drape Length: 6" / 15 cm	15" (38 cm)	100 / Box

BATRIK SHIELD™ FACE SHIELD

A cost-effective alternative to our premium SAVEFACE™ models, the Batrik Shield™ provides excellent splash protection. An ideal safety addition to all pandemic kits.

- Self-shaping foam head piece molds comfortably to forehead
- Shields can be used for longer periods due to comfort and improved materials
- Optical shield with anti-fog and anti-glare provide a clearer viewing area
- Soft plastic fasteners adhere the elastic band to the shield for additional comfort and easy use
- Latex free



PATENT PENDING

ASSEMBLED IN CANADA 

PRODUCT SELECTION


ITEM #	DESCRIPTION	DIMENSIONS	TOTAL LENGTH	PACKAGING
SF-096	BATRIK Shield™ no drape	Shield Length: 9" / 23 cm*	9" (23 cm)*	96 / Box
SF-896	BATRIK Shield™ with ANSI/AAMI PB70 Level 3 compliant protection drape	Shield Length: 9" / 23 cm* Drape Length: 7" / 17.8 cm	15.9" / 40.3 cm*	96 / Box

*Shield product tolerances: 8.9" - 9.5" / 22.5 cm - 23.2 cm

BATRIK SCRUB PAD

Our Scrub Pad is ideal for pre-cleaning your surgical instruments without leaving any metal residue that often occurs when using steel wool or metal sponges. The abrasiveness is perfectly balanced to provide you with efficient cleaning without damaging your instruments. Its fine fibers allow the user to polish, clean, and debur all types of material surfaces.

- Abrasive Material: Aluminum Oxide
- Grit: 400
- Grade: Very Fine
- Style: Flexible
- Pad material: Non-Woven
- Color: Maroon
- Type: Hand Pad



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
60-46050	Scrub Pad	4" x 6" / 10 cm x 15 cm	150 / Box (50 bags of 3 units each)
60-46003	Scrub Pad	4" x 6" / 10 cm x 15 cm	3 / Bag

Your First ™

PrAmLi Power³ Enzymes™ ENZYMATIC DETERGENTS

All PrAmLi Power³ Enzymes™ detergents are multi-enzyme pre-soak cleaners for medical equipment and instrumentation. Our proprietary unique blend of enzymes is effective in the removal of all proteins, blood, fatty tissue, and more. It will also remove stains from surgical scrubs, linens, and medical clothing.



Formulated with PrAmLi Power³ Enzymes™, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



PRODUCT SELECTION

ITEM #	DILUTION DETERGENT / WATER	PACKAGING
LOW FOAM DETERGENT		
SS-21-1001	4 to 8 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-1004	4 to 8 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
MILD FOAM DETERGENT		
SS-21-2001	4 to 8 ml / L (0.5 to 1oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-2004	4 to 8 ml / L (0.5 to 1oz / Gallon)	4 x 3.8 L (4 Gallons)
SUPER CONCENTRATE		
SS-21-3001	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	1 x 3.8 L (1 Gallon)
SS-21-3004	1 to 2 ml / L (0.125 to 0.25 oz / Gallon)	4 x 3.8 L (4 Gallons)

NOTE: Our Low Foam & Super Concentrate Detergent formulas are available in a 2.5 Gallon format based on a minimum order quantity. Contact us for more information.


GOOD TO KNOW:

OUR LOW FOAM DETERGENT MAY BE USED IN WASHERS/ ULTRASONIC CLEANERS. WE RECOMMEND CONFIRMING SUGGESTED DILUTION RATIO AND WATER TEMPERATURE WITH THE MANUFACTURER OF YOUR WASHER/ULTRASONIC CLEANER AS MACHINES/ SETTINGS VARY.

Your First ™

ENVIRO-SACHET™

ENZYMATIC GEL & NON-ENZYMATIC CLEANER

Your First ™ Enviro-Sachet™, our unique cleaning solutions now available in Enzymatic & Non-Enzymatic formats.



- Our Enzymatic Gel format is formulated with the same triple action **PrAmLi Power³ Enzymes™** that you know and trust.
- Our Non-Enzymatic Cleaner format is formulated with a free-rinsing, neutral pH detergent blend.

Both formats are powerful solutions that break down all bioburden for effective surface and instrument pre-cleaning.

PRODUCT SELECTION

ENZYMATIC GEL	NON-ENZYMATIC CLEANER	DESCRIPTION	PACKAGING
17-166	717-166	1 Dry Tubular Sponge & Enviro-Sachet™ set	25 / Box 4 Boxes / Case
17-266	717-266	1 Dry Flat Sponge & Enviro-Sachet™ set	50 / Box 4 Boxes / Case
17-177	717-177	1 Dry Tubular Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit
17-277	717-277	1 Dry Flat Sponge, Enviro-Sachet™ & Stainless Steel Bowl kit	1 Kit
17-EZSAMP	717-EZSAMP	SAMPLE KIT INCLUDES: 3 Dry Flat Sponges with Enviro-Sachet™ sets 1 Dry Tubular Sponge with Enviro-Sachet™ set 1 Reusable Stainless Steel Bowl	1 Kit
H-EZ500		Reusable Stainless Steel Bowl 4.3" x 5.3" x 3" / 11 cm x 13.5 cm x 7.5 cm	1 Unit

Your First ™

FOAMING SPRAY

READY-TO-USE ENZYME FOAMING SPRAY

Our Foaming Spray is ideal for instrument pre-cleaning as the thick foam spray enzymatically breaks down all bioburden to allow for effective cleaning and sterilization.

Formulated with **PrAmLi Power³ Enzymes™**, a unique blend of three powerful enzymes along with a surfactant to provide the ultimate in cleaning efficacy



INSTRUCTIONS



PRODUCT SELECTION

ITEM #	PACKAGING
SS-21-4001	1 x 946 ml (32 oz) Spray Bottle
SS-21-4012	12 x 946 ml (32 oz) Spray Bottles

Important: In addition to our enzymatic formulation, our spray bottle contains as much as 25% more solution when compared to the competition.

Your First ™

FIRST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-A0909	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case
17-A1212	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case

Your First ™

DRY SPONGES

Our sponges are specifically fabricated for use with all medical instrumentation. Unlike most competitor sponges, our sponges are white to better evidence the effectiveness of the cleaning process. Dark colored sponges easily hide contaminants, but with our white sponges you can visually evaluate the pre-clean efficacy.



Available in two formats, Flat and Tubular:

Will not flake or scratch expensive equipment

Eco-friendly

100% medical grade high-density foam

Compatible with enzymatic and disinfectant solutions

TIP:
OUR TUBULAR SPONGE CAN ALSO BE USED AS A BUMPER TO PROTECT YOUR SCOPES

High absorbency retains more solution than conventional sponges or paper wipes

OUR WHITE SPONGES MAY APPEAR SLIGHTLY YELLOWED. THIS IS ONLY DUE TO UV EXPOSURE, IS NOT A REFLECTION OF ITS AGE.

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-100	DRY TUBULAR SPONGE White Dispenser Box	25 / Box (4 Boxes / Case)
17-101	DRY TUBULAR SPONGE Bulk	560 / Case
17-200	DRY FLAT SPONGE White Dispenser Box	50 / Box (4 Boxes / Case)
17-201	DRY FLAT SPONGE Bulk	1,000 / Case

Your First ™

ENZYMATIC BEDSIDE SPONGE KITS



Our Enzymatic Bedside Sponge Kits feature sponges all pre-treated with **PrAmLi Power³ Enzymes™** to provide the most effective pre-clean or surface clean solution. This unique three-enzyme blend provides you with confident and effective cleaning. The tubular sponge is ideal for all endoscopes.




PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
17-111	DRY TUBULAR SPONGE Moisten in any 250-500 ml basin to activate the enzymatic concentrate.	25 / Box (4 Boxes / Case)
17-211	DRY FLAT SPONGE Moisten in any 250-500 ml basin to activate the enzymatic concentrate.	50 / Box (4 Boxes / Case)
17-122	DRY TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	15 / Box (4 Boxes / Case)
17-222	DRY FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	25 / Box (4 Boxes / Case)
17-133	READY TO USE TUBULAR SPONGE Comes bathing in 500 ml of enzymatic solution.	24 / Case
17-233	READY TO USE FLAT SPONGE Comes bathing in 500 ml of enzymatic solution.	24 / Case
17-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of enzymatic solution.	12 / Box (4 Boxes / Case)
17-155	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case
17-255	READY TO USE FLAT SPONGE Comes bathing in 220 ml of enzymatic solution.	24 / Case

Your First ™

NON-ENZYMATIC DETERGENT BEDSIDE SPONGE KITS

Our Bedside Sponge Kits with Non-Enzymatic Detergent are **Your First **™ in the medical instrument cleaning process.



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
717-122	DRY TUBULAR SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	15 / Box (4 Boxes / Case)
717-222	DRY FLAT SPONGE Comes in a freestanding pouch to use as a basin. Simply add 500 ml of water.	25 / Box (4 Boxes / Case)
717-133	READY TO USE TUBULAR SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
717-233	READY TO USE FLAT SPONGE Comes bathing in 500 ml of detergent solution.	24 / Case
717-144	READY TO USE DUAL PURPOSE TUBULAR SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
717-244	READY TO USE DUAL PURPOSE FLAT SPONGE Comes in a freestanding pouch soaked in 100 ml of detergent solution.	12 / Box (4 Boxes / Case)
717-155	READY TO USE TUBULAR SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case
717-255	READY TO USE FLAT SPONGE Comes bathing in 220 ml of detergent solution.	24 / Case

Your First ™

TRANSPORT BAGS FOR ENDOSCOPE AND MEDICAL INSTRUMENTS

Thoughtfully designed to eliminate touch-points during endoscope transport from point-of-care to the reprocessing area, and back, focused on safety of patients and medical staff.

- Made of recyclable LDPE material known to be quite flexible, highly resilient, and tough
- Latex and DEHP-free
- Excellent resistance to dilute and concentrated acids, alcohols, bases and esters*
- Good resistance to aldehydes, ketones, and vegetable-based oils*
- Can withstand high temperatures of up to 90 °C (194°F)
- Tamper-evident seal for added safety, also aids in keeping the instrument moist to avoid drying of bio-debris
- Clearly marked bags with color-coded visuals facilitates immediate identification of clean vs. used instruments
- Users can write important information directly on the bag in spaces marked with visual icons
- Available in two convenient sizes

GOOD TO KNOW: USING **Your First ™** SEALABLE TRANSPORT BAGS REDUCES USER ERROR, PREVENTS CROSS-CONTAMINATION, AND PROMOTES PATIENT AND MEDICAL STAFF SAFETY.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING
17-HAP24	Clean Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
17-SAD24	Soiled Instruments Bag	24" x 22" / 61.0 cm x 56.0 cm	200 / Case
17-HAP30	Clean Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case
17-SAD30	Soiled Instruments Bag	30.5" x 26" / 77.0 cm x 66.0 cm	200 / Case

*Reference on file

CLEAN INSTRUMENTS BAG



SOILED INSTRUMENTS BAG



MADE IN CANADA

BATRIK ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.

In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

RIGID



FLEXIBLE



12" / 30 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-32504-S	411-32504	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-32505-S	411-32505	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-35005-S	41-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35006-S	41-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35008-S	41-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35010-S	41-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35012-S	41-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35015-S	41-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-35020-S	41-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-46502-S	41-46502	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-47503-S	41-47503	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-47504-S	411-47504	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-47505-S	411-47505	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-50005-S	41-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50006-S	41-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50008-S	41-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50010-S	41-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50012-S	41-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50015-S	41-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-50020-S	41-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH



	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-61502-S	41-61502	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	41-62503-S	41-62503	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	411-62504-S	411-62504	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	411-62505-S	411-62505	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
RIGID BRUSH & SHEATH	41-65005-S	41-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	41-65006-S	41-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	41-65008-S	41-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65010-S	41-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65012-S	41-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65015-S	41-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-65020-S	41-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

30" / 75 CM SINGLE-END WIRE BRUSH


	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-76502-S	41-76502	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-77503-S	41-77503	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-77504-S	411-77504	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-77505-S	411-77505	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-80005-S	41-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80006-S	41-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80008-S	41-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80010-S	41-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80012-S	41-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80015-S	41-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-80020-S	41-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3



FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm



RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

2mm 0.078"	3 mm 0.118"	4 mm 0.157"	5 mm 0.196"	6 mm 0.236"	8 mm 0.314"	10 mm 0.394"	12 mm 0.472"	15 mm 0.591"	20 mm 0.787"

BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES

ASYMMETRIC

14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES SYMMETRIC

14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

ASYMMETRIC BRUSHES

5 MM
6 MM

8 MM
10 MM

10 MM
12.5 MM

15 MM
20 MM



BATRIK

REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED WIRE BRUSHES FEATURE A PROTECTIVE BLUE TIP FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
21-10039	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
21-10067	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
20-10059	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
20-10079	6 FR	1" / 25.4 mm	0.079" / 2 mm
20-10090	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	1" / 25.4 mm	0.118" / 3 mm
20-10129	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
20-10146	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
20-10157	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3

12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
22-12078	2" / 50.8 mm	0.078" / 2 mm
22-12118	2" / 50.8 mm	0.118" / 3 mm
22-12157	3" / 76 mm	0.157" / 4 mm
22-12197	3" / 76 mm	0.197" / 5 mm
22-12236	3" / 76 mm	0.236" / 6 mm
22-12393	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3

16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
23-16197	2.5" / 63.5 mm	0.197" / 5 mm
23-16236	2.5" / 63.5 mm	0.236" / 6 mm
23-16276	3" / 76 mm	0.276" / 7 mm
23-16394	3" / 76 mm	0.394" / 10 mm
23-16472	3" / 76 mm	0.472" / 12 mm
23-16591	3" / 76 mm	0.591" / 15 mm
23-16787	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3

18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
19-18010	0.47" / 12 mm	0.039" / 1 mm
19-18015	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3

24" / 610 MM NYLON TWISTED WIRE BRUSHES


ITEM #	BRUSH LENGTH	BRUSH (Ø)
24-24098	2.5" / 63.5 mm	0.098" / 2.5 mm
24-24197	2.5" / 63.5 mm	0.197" / 5 mm
24-24276	3" / 76 mm	0.276" / 7 mm
24-24394	3" / 76 mm	0.394" / 10 mm
24-24472	3" / 76 mm	0.472" / 12 mm

PACKAGING: Bag of 3

TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
25-08500	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3



BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3

BATRIK STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH
TECHNIQUE IS HIGHLY
RECOMMENDED WITH THESE
BRUSHES. SIMPLY INSERT THE
BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

BATRIK INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



50-40919



50-40922



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



50-40920



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



50-40918



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40917



50-40921


PRODUCT SELECTION

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20

BATRIK REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT DESCRIPTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	3 / Bag

BATRIK

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES THE NEED FOR THE STAINLESS STEEL TOOTHBRUSH-STYLE BRUSH, WHICH CAN DAMAGE VALUABLE INSTRUMENTS!

PRODUCT SELECTION

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm

PACKAGING: Bag of 3

BATRIK REUSABLE INSTRUMENT CLEANING BRUSHES


Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH



STERILIZER DRAIN BRUSH



PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm

PACKAGING: Bag of 3

BATRIK ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72" / 120 mm Ø 0.47" / 12 mm	1.18" / 30 mm Ø 0.3125" / 8 mm	2" / 50 mm Ø 1.75" / 45 mm

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

BATRIK ALL PURPOSE REUSABLE BRUSHES

Plastic Handle Brushes

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.



28-72550



28-69394



29-72550



29-69394



29-70050



50-40911

PRODUCT SELECTION

ITEM #	DESCRIPTION	LENGTH TOTAL	BRUSH LENGTH	BRUSH HEIGHT
28-72550	Plastic Single End Nylon Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-72550	Plastic Single End Stainless Steel Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-70050	Plastic Single End Brass Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
28-69394	Plastic Double End Nylon Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
29-69394	Plastic Double End Stainless Steel Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
50-40911	Hand Scrub Nylon Brush Reusable autoclavable	4.09" / 104 mm	3.74" / 95 mm	1.1" / 28 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 8.89 cm x 20.32 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 11.5 cm x 20.32cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm
		SF-100	2 SAVEFACE™ Face Shields	9" / 23 cm

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-100

Packaging: 100 / Box

BATRIK HYSTEROSONOGRAPHY CATHETER

Hystero-sonography (also known as sonohysterography) can now be performed safely, simply and atraumatically with the gentle infusion of saline into the uterus with the use of our H/S Catheter.


- Available in two sizes: 5FR and 7FR
- Minimally invasive and ideal for office use
- Gentle placement sheath included (decreasing the need for a tenaculum) to minimize patient discomfort
- Flexible and easy-to-use
- Total length of 410mm with a working length of 280 mm
- Latex Free

GOOD TO KNOW: ACHIEVE
IMPROVED IMAGERY
AND DIAGNOSIS WITH
SONOHYSTEROGRAMS

“ We depend upon the reliability, precision and functionality of the Batrik H/S Catheter in our clinical and research operations.

Having performed more than 100,000 sonohysterography procedures, we can unequivocally say that this is the best HSG catheter that we have ever used

”
-Dr. Alex Hartman



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
51-5005	5 FR H/S Catheter	10 / box
51-5007	7 FR H/S Catheter	10 / box

SURGMED

GROUP ● ● ●

Proudly Canadian, **The Surgmed™ Group of Companies** began operations in 1992 as a single distribution company known as **Advanced Surgi-Pharm Inc.** The company quickly became known as an industry leader and specialty distributor of value-added surgical products in the Canadian hospital market. Five years later, **Surgmed expanded upon its success with the launch of Batrik Medical Manufacturing Inc.** in 1997 and the subsequent acquisition of **Imperial Surgical Ltd.** one year later in 1998.


Together, our three unique companies serve the Canadian and global healthcare markets with high quality and competitively priced products, all backed by world-class customer care.



Founded in 1992, Advanced Surgi-Pharm Inc. is the distribution arm of Surgmed™ in Canada, representing leading US, international and domestic manufacturers for a wide range of specialty products with unique benefits for applications in the Operating Room, Women’s Health, Labor & Delivery, ICU, and Sterile Processing areas. Equipped with a direct sales team, ASP provides effective nationwide coverage as the Canadian distributor of Batrik, Vygon, Rocket Medical, Cooper Surgical, SunMed, and other industry leading products including a wide-ranging portfolio of top-rated PPE products.

Established in 1997 and sold in 80 countries globally, Batrik Medical Manufacturing Inc. specializes in the manufacturing of endoscopy accessories, enzymatic solutions, infection control, central sterilization, and surgical devices. Batrik’s in-house controlled environment production facility adheres to strictly enforced quality and regulatory control systems. They are registered as a device manufacturing establishment with the United States FDA and are in full compliance with CE and Health Canada Regulations. All Batrik products are rigorously tested, liability insured, and meet/exceed high international quality standards.

Established in 1935, Imperial Surgical Ltd. is one of Canada’s oldest hospital capital equipment manufacturers, well-known for their stainless steel healthcare equipment for Sterile Processing, Operating Room, Examination and Treatment, Patient Handling, and Mortuary/ Pathology areas. Their growing portfolio of products also has applications in Infection Control & Prevention, an important area of focus for their future.



SURGMED GROUP

850 Halpern Ave., (Dorval) Montreal QC H9P 1G6
Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax.: +1 514-631-9083
Email: info@surgmed.com | Website: www.batrik.com