

# BOWL STANDS & RECEPTACLES

Our Bowl Stands and Receptacles are designed for convenience, mobility, and durability. Available in five versatile models that accommodate most standard sizes of stainless steel bowls and pails. Fully welded stainless steel construction throughout with a #4 polished finish.

All bowl stands are backed by an outstanding **7 year warranty** on all parts (see warranty for details).

## KEY FEATURES

- Seamless welded tubular frame with lower cross-bracing provides maximum rigidity and extended product life
- Optional lower stainless steel shelf is available for all configurations
- Non-marking conductive swivel design 2" (51 mm) casters

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>OR-3631</b>	<b>SPONGE RECEPTACLES: FLOOR UNITS</b> Mobile Kick Bucket with Qt. Pail and Perimeter Bumper Receptacle: 14.5" Ø x 8" H / 368 mm Ø x 203 mm H Receptacle and bucket: 14.5" Ø x 15" H / 368 mm Ø x 381 mm H Bucket: 12" Ø x 11" H / 305 mm Ø x 279 mm H
<b>OR-3632</b>	<b>MOBILE SPONGE RECEPTACLES</b> with 10 Qt. Oval Basin and Perimeter Bumper
<b>OR-3631-TB</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Mobile Table Height Receptacle with 20 Qt. Pail 15" L x 16" W x 38" H (381 mm x 406 mm x 965 mm)
<b>OR-3651</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Single Bowl Stand with 7 Qt. Basin 15" L x 16" W x 34" H (381 mm x 406 mm x 864 mm)
<b>OR-3652</b>	<b>SPONGE RECEPTACLES: TABLE HEIGHT</b> Double Bowl Stand with Two 7 Qt. Basins 33" L x 16" W x 34" H (838 mm x 406 mm x 864 mm)



THE LOWER SHELF IS AN OPTIONAL FEATURE THAT MUST BE ADDED TO THE ITEM IF DESIRED

## ACCESSORIES

ITEM #	DESCRIPTION
<b>OR-3631-PS</b>	12 QT. STAINLESS STEEL PAIL for OR-3631
<b>OR-3631-TB-PS</b>	20 QT. STAINLESS STEEL PAIL for OR-3631-TB
<b>OR-3632-BS</b>	10 QT. OVAL STAINLESS STEEL BASIN for OR 3632
<b>OR-3651-BS</b>	7 QT. STAINLESS STEEL BASIN for OR-3651 and OR-3652

Add suffix -LS to any Item # to add a lower shelf