

CATHETER & INSTRUMENT TRAYS



H-0153




H-0156



H-0159

PRODUCT SELECTION

ITEM #	DESCRIPTION
H-0151	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0152	COVER FOR H-0151
H-0153	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 219 mm L x 81 mm W x 35 mm H
H-0154	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0155	COVER FOR H-0154
H-0156	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 231 mm L x 132 mm W x 49 mm H
H-0157	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0158	COVER FOR H-0157
H-0159	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 320 mm L x 95 mm W x 50 mm H



CATHETER & INSTRUMENT TRAYS

PRODUCT SELECTION

ITEM #	DESCRIPTION
H-0160	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0160P	PERFORATED CATHETER & INSTRUMENT TRAY (NO COVER)
H-0161	COVER FOR H-0160 AND H-0160P
H-0162	CATHETER & INSTRUMENT TRAY SET (WITH COVER, FLUSH HANDLE) 310 mm L x 196 mm W x 50 mm H
H-0163	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0163P	PERFORATED CATHETER & INSTRUMENT TRAY (NO COVER)
H-0164	COVER FOR H-0163 AND H-0163P
H-0165	CATHETER & INSTRUMENT TRAY SET (WITH COVER, FLUSH HANDLE) 310 mm L x 196 mm W x 65 mm H
H-0166	CATHETER & INSTRUMENT TRAY (NO COVER)
H-0167	COVER FOR H-0166
H-0168	CATHETER & INSTRUMENT TRAY SET (WITH COVER, RAISED HANDLE) 441 mm L x 116 mm W x 45 mm H



H-0162



H-0165



H-0168

