

# CHAIRS & STOOLS

Our Chairs & Stools revolve and are an economical solution for an operating room or clinical application. Available with and without upholstery in four models. Stainless steel construction facilitates easy clean-up and extended product life. Fully welded stainless steel seamless tubular frame construction with a #4 polished finish.

All chairs & stools are backed by our **1 year warranty** on structural and all components (excluding casters).

## KEY FEATURES

- Wide base provides added stability especially when used with optional casters
- Stainless steel 13" (330 mm) diameter seat with optional backrest
- Upholstered 14" (356 mm) diameter seats feature 2" (51 mm) dual density foam core and antibacterial conductive double stitched vinyl cover
- Integrally welded stainless steel foot ring for maximum strength and comfort
- Heavy-duty threaded spindle and bracket assembly provide easy height adjustment and rigid support
- Choose between non-marking conductive swivel design 2" (51 mm) casters or conductive rubber crutch tips



OR-3620



OR-3622

## PRODUCT SELECTION

| ITEM #         | DESCRIPTION                                                                                                        |
|----------------|--------------------------------------------------------------------------------------------------------------------|
| <b>OR-3620</b> | STAINLESS STEEL REVOLVING STOOL<br>Height adjustment 19" to 27" (483 mm to 686 mm)                                 |
| <b>OR-3621</b> | STAINLESS STEEL REVOLVING STOOL<br>Height adjustment 25" to 32" (635 mm to 813 mm)                                 |
| <b>OR-3622</b> | REVOLVING STOOL WITH UPHOLSTERED SEAT, AND STAINLESS STEEL BASE<br>Height adjustment 19" to 27" (483 mm to 686 mm) |
| <b>OR-3623</b> | REVOLVING STOOL WITH UPHOLSTERED SEAT, AND STAINLESS STEEL BASE<br>Height adjustment 25" to 32" (635 mm to 813 mm) |

## ACCESSORIES

| ITEM #         | DESCRIPTION                                                                          |
|----------------|--------------------------------------------------------------------------------------|
| <b>OR-3628</b> | CASTER ASSEMBLY WITH SOCKET KIT                                                      |
| <b>OR-3629</b> | CONDUCTIVE SEAT COVER, AND PAD FOR 13" / 330 MM DIAMETER STAINLESS STEEL STOOL SEATS |