

# CLOSED CASE CARTS



STREAMLINE YOUR WORKFLOW  
WITH OUR CLOSED CASE CARTS  
**PROUDLY MADE IN CANADA** 



# CLOSED CASE CARTS

Our Closed Case Carts come in a wide range of sizes with flexible shelving options, ergonomic doors, smart design choices, safety features, superior maneuverability, and optional accessories. They are robust and durable yet lightweight thanks to our one-piece body-wrap design of 18-gauge stainless steel.

All closed case carts are backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

## KEY FEATURES

### SMART DESIGN

- Top surface has a drip-proof ledge for spill drainage, while 5 degree sloped bottom pan and perforated chassis aid in drainage.
- The external welded handle maximizes interior storage.
- The outer door lock feature allows you to keep an open door position during washing, loading and unloading.



### ERGONOMIC DOOR ASSEMBLIES

- Strong double pan design features bottom edge drainage holes (outside positioning)
- Heavy gauge hinge assemblies are easily removed for cleaning/adjustment
- SteriLock spring loaded door latch prevents accidental opening during transportation
- Include external push-pull handle for safe locking



### OPTIONAL FEATURES

- Choose optional card holders, paddle latch assembly and retractable tow system assembly.

### FLEXIBLE SHELVING

- Adjustable open shelf design provides additional shelf width and depth for large containers.
- Shelves types include solid, perforated or wired slide out featuring roll-out integral stops to prevent tipping during use.

### SAFETY FEATURES

- Full perimeter bumper made from high temperature resistant rubber provides impact protection.
- Permanently sealed seams and crevices prevents water accumulation and sources of contamination.

## SUPERIOR MANEUVERABILITY FEATURES

- Each cart features two all-swivel, one steer-lock, and one total-lock caster to facilitate easy steering, handling, and positioning.
- Stainless steel caster body is cart washer compatible and provides a high load rating of 225 lbs per wheel.



## ACCESSORIES

| ITEM #           | DESCRIPTION   |
|------------------|---|
| <b>OR-3875</b>   | CARD HOLDER<br>6" W x 4"H (152mm x 102mm)                                 |
| <b>OR-3876</b>   | CARD HOLDER<br>9" W x 11" H (229mm x 279mm)                               |
| <b>OR-3877</b>   | SECURITY / RING-SEAL BRACKETS   |
| <b>OR-3881</b>   | PUSH / PULL HANDLE  |
| <b>OR-3871-S</b> | PADDLE LATCH ASSEMBLY SMALL   |
| <b>OR-3871-L</b> | PADDLE LATCH ASSEMBLY LARGE   |
| <b>OR-3840</b>   | SOLID SLIDE-OUT SHELF<br>24" W x 24" H - Only for single door models      |
| <b>OR-3841</b>   | PERFORATED SLIDE-OUT SHELF<br>24" W x 24" H - Only for single door models |
| <b>OR-3843</b>   | WIRE SLIDE-OUT SHELF<br>24" W x 24" H - Only for single door models       |
| <b>OR-3845</b>   | SOLID SLIDE-OUT SHELF<br>24" W x 36" H - Only for double door models      |
| <b>OR-3846</b>   | PERFORATED SLIDE-OUT SHELF<br>24" W x 36" H - Only for single door models |
| <b>OR-3848</b>   | WIRE SLIDE-OUT SHELF<br>24" W x 36" H - Only for single door models       |



## PRODUCT SELECTION

| ITEM #         | DESCRIPTION                      | DIMENSIONS<br>(L x W x H) IN / MM   | CAPACITY<br>(CU. FT.) * |
|----------------|----------------------------------|---|-------------------------|
| <b>OR-3810</b> | SINGLE DOOR /<br>STANDARD HEIGHT | Shelf: 24" x 24" (610 mm x 610 mm)<br>Cabinet: 27" x 25 1/2" x 33" (686 mm x 648 mm x 838 mm)<br>Overall: 33 1/4" x 28" x 39 1/2" (845 mm x 711 mm x 1003 mm)   | 11.5                    |
| <b>OR-3815</b> | SINGLE DOOR /<br>TALL DESIGN     | Shelf: 24" x 24" (610 mm x 610 mm)<br>Cabinet: 27" x 25 1/2" x 48" (686 mm x 648 mm x 1219 mm)<br>Overall: 33 1/4" x 28" x 54 1/2" (845 mm x 711 mm x 1384 mm)  | 16.5                    |
| <b>OR-3820</b> | DOUBLE DOOR /<br>STANDARD HEIGHT | Shelf: 36" x 24" (914 mm x 610 mm)<br>Cabinet: 39" x 25 1/2" x 33" (991 mm x 648 mm x 838 mm)<br>Overall: 45 1/4" x 28" x 39 1/2" (1149 mm x 711 mm x 1003 mm)  | 16.6                    |
| <b>OR-3825</b> | DOUBLE DOOR / TALL<br>DESIGN     | Shelf: 36" x 24" (914 mm x 610 mm)<br>Cabinet: 39" x 25 1/2" x 48" (991 mm x 648 mm x 1219 mm)<br>Overall: 45 1/4" x 28" x 54 1/2" (1149 mm x 711 mm x 1384 mm) | 24.2                    |



**SURGMED GROUP**

2165 46<sup>th</sup> Ave., Lachine, Quebec H8T 2P1

Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax: +1 514-631-9083

Email: [info@surgmed.com](mailto:info@surgmed.com) | Website: [www.imperialsurgical.com](http://www.imperialsurgical.com)