

CLOSED CASE CARTS

Our Closed Case Carts are built to last and come in a wide range of sizes with flexible shelving options, ergonomic door assemblies, smart design choices, important safety features, superior maneuverability, and optional accessories. They are highly robust and durable yet lightweight thanks to our one-piece body-wrap design of 18-gauge superior stainless steel.

All closed case carts are backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H) IN / MM	CAPACITY (CU. FT.) *
OR-3810	SINGLE DOOR / STANDARD HEIGHT	Shelf: 24" x 24" (610 mm x 610 mm) Cabinet: 27" x 25 1/2" x 33" (686 mm x 648 mm x 838 mm) Overall: 33 1/4" x 28" x 39 1/2" (845 mm x 711 mm x 1003 mm)	11.5
OR-3815	SINGLE DOOR / TALL DESIGN	Shelf: 24" x 24" (610 mm x 610 mm) Cabinet: 27" x 25 1/2" x 48" (686 mm x 648 mm x 1219 mm) Overall: 33 1/4" x 28" x 54 1/2" (845 mm x 711 mm x 1384 mm)	16.5
OR-3820	DOUBLE DOOR / STANDARD HEIGHT	Shelf: 36" x 24" (914 mm x 610 mm) Cabinet: 39" x 25 1/2" x 33" (991 mm x 648 mm x 838 mm) Overall: 45 1/4" x 28" x 39 1/2" (1149 mm x 711 mm x 1003 mm)	16.6
OR-3825	DOUBLE DOOR / TALL DESIGN	Shelf: 36" x 24" (914 mm x 610 mm) Cabinet: 39" x 25 1/2" x 48" (991 mm x 648 mm x 1219 mm) Overall: 45 1/4" x 28" x 54 1/2" (1149 mm x 711 mm x 1384 mm)	24.2

ACCESSORIES

ITEM #	DESCRIPTION
OR-3875	CARD HOLDER 6" W x 4"H (152mm x 102mm)
OR-3876	CARD HOLDER 9" W x 11" H (229mm x 279mm)
OR-3877	SECURITY / RING-SEAL BRACKETS
OR-3881	PUSH / PULL HANDLE
OR-3871-S	PADDLE LATCH ASSEMBLY SMALL
OR-3871-L	PADDLE LATCH ASSEMBLY LARGE
OR-3840	SOLID SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3841	PERFORATED SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3843	WIRE SLIDE-OUT SHELF 24" W x 24" H - Only for single door models
OR-3845	SOLID SLIDE-OUT SHELF 24" W x 36" H - Only for double door models
OR-3846	PERFORATED SLIDE-OUT SHELF 24" W x 36" H - Only for single door models
OR-3848	WIRE SLIDE-OUT SHELF 24" W x 36" H - Only for single door models



OR-3825

KEY FEATURES

FLEXIBLE SHELVING

- Adjustable open shelf design provides additional shelf width and depth for large containers
- Shelves types include solid, perforated, or wire with optional roll-out featuring integral stops to prevent tipping during use

ERGONOMIC DOOR ASSEMBLIES

- Strong double pan design features bottom edge drainage holes (outside positioning)
- Heavy gauge hinge assemblies are easily removed for cleaning/adjustment
- SteriLock spring loaded door latch prevents accidental opening during transportation and ensures a fully open position during washing, loading, and unloading



SAFETY FEATURES

- Full perimeter bumper made from high temperature resistant rubber provides impact protection
- Permanently sealed seams and crevices prevents water accumulation and sources of contamination

SMART DESIGN

- Top surface features an integral drip-proof ledge to ensure proper drainage of accidental spills without contaminating cart interior
- Cabinet bottom pan is sloped 5 degrees and chassis support frame is perforated to facilitate drainage
- Externally mounted integrally welded push-pull handle maximizes interior storage space

SUPERIOR MANEUVERABILITY

- Each cart features two all-swivel, one steer-lock, and one total-lock caster to facilitate easy steering, handling, and positioning
- Stainless steel caster body is cart washer compatible and provides a high load rating of 225 lbs per wheel
- Optional retractable linkage/training assembly provides for easy manual cart training/towing

OPTIONAL ACCESSORIES

- Choose amongst optional card holders, paddle latch assembly, linkage, and roll-out shelving

