

PROTECTIVE LAMINATED STITCHED GOWNS

NON-STERILE
STITCHED SEAMS

Our Disposable Protective Laminated Gowns are impermeable to fluids, offering exceptional protection to healthcare workers.

- Fabric Compliant to EN 20811, EN ISO 22610 and ASTM F739-12
- 56 gr/m² fabric for maximum **protection, breathability, and comfort**
- Knitted soft wrist cuffs
- Tie closure
- **Stitched seams**
- **Non-Sterile**



PRODUCT SELECTION

ITEM #	SIZE	OVERALL LENGTH	SLEEVE LENGTH	DIAMETER
SF-302-S	S	110 cm / 43"	55 cm / 22"	140 cm / 55"
SF-302-M	M	115 cm / 45"	58 cm / 23"	145 cm / 57"
SF-302-L	L	127 cm / 50"	60 cm / 24"	150 cm / 59"
SF-302-XL	XL	135 cm / 53"	62 cm / 24"	150 cm / 59"
SF-302-XXL	XXL	140 cm / 55"	64 cm / 25"	155 cm / 61"
SF-302-XXXL	XXXL	145 cm / 57"	66 cm / 26"	160 cm / 63"

PACKAGING: 50 GOWNS / BOX