

EQUIPMENT TABLES

Our Equipment Tables come in three sizes and are a sterile and efficient solution for organizing and storing equipment and supplies. Fully welded heavy gauge stainless steel construction is easy to clean and maintain.

All equipment tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- All models include stainless steel welded drawer liners, a lower shelf, and three -sided ribbon rail
- Heavy-duty conductive swivel casters come in 2" (51 mm) or 3" (76 mm) depending on table size



PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
OR-3644	EQUIPMENT/INSTRUMENT TABLE WITH TWO DRAWERS, GUARD RAIL, AND LOWER SHELF 36" x 18" x 34" (914 mm x 457 mm x 864 mm)
OR-3648	EQUIPMENT/INSTRUMENT TABLE WITH ONE DRAWER, GUARD RAIL, AND LOWER SHELF 20" x 18" x 34" (508 mm x 457 mm x 864 mm)
OR-3650	ANESTHETIST'S TABLE WITH THREE DRAWERS (2 X 4"/102MM DEEP, 1 X 8"/204MM DEEP), GUARD RAIL, AND LOWER SHELF 20" x 18" x 34" (508 mm x 457 mm x 864 mm)