

# FIRST AID COUCH

Our First Aid Couch is a general-purpose treatment couch with a full-length vinyl upholstered top section featuring a contoured head support area and a dual density foam core for maximum comfort and support.

## KEY FEATURES

- Vinyl material is sulphide stain resistant and meets standard flame resistance requirements
- Full welded tubular leg frame assembly with height adjustable levelling feet (sub-assembly required)
- Custom colours available: Aqua Green, Colonial Blue, Mauve (note: all custom colours may require a set-up charge, subject to availability and volume)



ET-1574

# EXAMINATION CHAIRS & STOOLS

Our general-purpose examination chairs and stools feature a pneumatic assist gas cylinder and stable five leg base.

## KEY FEATURES

### BASE

- Five leg, 22"/ 559 mm diameter base assembly
- Reinforced nylon resin base with steel collar for maximum stability

### CASTERS

- 2"/51 mm dual wheel general purpose casters for tile or carpet surfaces

### SEAT HEIGHT

- Adjustable by pneumatic cylinder with dual nylon wheel

### UPHOLSTERY

- Double density 3"/76 mm foam upholstered top
- 15 ¾"/ 400 mm diameter seat
- Vinyl material is sulphide stain resistant with a high abrasion resistance factor and meets standard flame resistance requirements



ET-1955

## OPTIONAL FEATURES

- Three cylinder seat height configurations from 16"/406 mm to 32"/ 813 mm
- Adjustable backrest depth is fully adjustable (moved front or back)
- Adjustable footring assembly for tall cylinder configurations
- Safety Casters

## PRODUCT SELECTION

ITEM #	DESCRIPTION
ET-1574	FIRST AID COUCH
ET-1955	STANDARD HEIGHT EXAM CHAIR/STOOL
ET-1960	TALL EXAM CHAIR/STOOL
ET-1961	ADJUSTABLE BACKREST FOR STOOL
ET-1970	LOCKING CASTERS set of 5

All couches, chairs & stools are backed by our **1 year warranty** on structural and all components (excluding casters).