

FOOT STOOLS

Our Foot Stools are available in a variety of configurations for a full range of applications in the operating room. Durable stainless steel fully welded construction provides maximum strength and facilitates easy clean-up and extended product life.

All foot stools are backed by our **1 year warranty** on structural and all components.

KEY FEATURES

- Single or double non-slip rubber tread surfaces with flared leg design for superior stability
- Heavy gauge stainless steel platform base with die-formed edge
- Seamless welded tubular frame leg support with conductive rubber crutch tips
- Fully welded assembly



PRODUCT SELECTION

ITEM #	DIMENSIONS (L x W x H)
OR-36360	16" x 12" x 8" (406 mm x 305 mm x 203 mm)
OR-36360-B	24" x 12" x 8" (610 mm x 305 mm x 203 mm)
OR-36360-E	30" x 12" x 8" (762 mm x 305 mm x 203 mm)

Add suffix -H to Item # to add a Stainless Handrail