


LINEN INSPECTION TABLES

Our Linen Inspection Tables provide a large, LED illuminated work surface for the inspection of linen drapes and packing wraps. Available in 2 sizes and constructed from stainless steel with a sealed translucent plexiglass work surface. The stainless steel frame and legs ensures durability for extended product life.

All linen inspection tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- LED fixtures provide bright, consistent lighting, are power efficient and last longer than traditional bulbs
- Optional lower shelf provides extra storage space
- Tubular legs are available with optional heavy-duty grey non-marking casters with brakes
- Features height adjustable levelers
- CSA Approved hospital grade power cord (10 feet / 3 meters) and light switch
- Easily accessible lamp assembly for bulb replacement
- Service required: 120V, 60Hz



PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS (L x W x H)
SP-8210 LS	LINEN INSPECTION TABLE WITH LOWER SHELF (Standard Length)	72"x 32"x 36" (1829 mm x 813 mm x 914 mm)
SP-8215	CASTER KIT FOR STANDARD TABLE (Set of Four) Factory Installed	
SP-8220 LS	LINEN INSPECTION TABLE WITH LOWER SHELF (Standard Length)	
SP-8225	CASTER KIT FOR EXTENDED TABLE (Set of Six)	