



# STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, our patented Stackable Step Stools have been designed with infection prevention as a prime focus.



## PRODUCT FEATURES

- Weight: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Transport Cart Available
- Weight Capacity: 300 lb (136 kg)
- Warranty: 7 years (step stool only)

## ANTIMICROBIAL COATING:

Prevents the spread of hazardous microorganisms and suppresses the growth of mold, bacteria, yeast, and mildew.

## ANTI-SLIP SURFACE:

Provides enhanced stability.

## RUBBER-FREE CORNERS:

Avoids fluid collection.

## WASHABLE CART:

Model OR-36364 (without mat) can be sterilized.

## HYGIENIC DESIGN:

No welds or large hollow openings to trap debris or fluids.

## PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR
OR-36363	Stackable Step Stool with mat	Black
OR-36363-02	Stackable Step Stool with mat	Dark Blue
OR-36363-03	Stackable Step Stool with mat	Green
OR-36363-07	Stackable Step Stool with mat	Orange
OR-36364	Stackable Step Stool without mat	N/A

Packaging: Sold by Unit

## REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING
OR-36363-1A	Replacement Mat for OR-36363	1 Unit
OR-36363-1B	Replacement Foot Pegs	1 Unit

## OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING
OR-3693	Transport cart	1 Unit
809-1705	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box
899-1710	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box



OR-3693