

At Imperial Surgical Ltd., we have two types of IV Poles to offer. Our versatile modular design means you simply select the type and the components that work best for your needs. All of our IV Poles and components are backed by an industry leading 10 year warranty against manufacture defects.

CHOOSE AMONGST

- Our traditional **Stainless Steel IV Poles**, made from 300 grade stainless steel that is easy to clean and disinfect
- Our **NEW Antimicrobial IV Poles**, made from the same high-quality stainless steel and protected with PMRA and EPA registered **CuVerro Shield™** by **Aereus Technologies**

STAINLESS STEEL IV POLES

ITEM #	DESCRIPTION
OR-36770	IV POLE, CENTRE COLLAR, AND KNOB with 23" cream aluminium base (V-400016), and 2" casters
OR-36760-25	IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-400013-25), and 2" casters
OR-36760-253	with 25" stainless steel base (V-400013-25), and 3" casters
OR-36771-25	IV POLE, CENTRE COLLAR, AND KNOB with 25" stainless steel base (V-40013-25), and 2" casters, and weight kit of 13.5 lbs. (V-50004-25)
OR-36771-253	with a 25" stainless steel base (V-40013-25), and 3" casters, and weight kit of 13.5 lbs. (V-50004-25)
OR-36772-253-UHN	IV POLE, CENTRE COLLAR, AND KNOB with 25" steel base (V-400014-25-ST), and 3" casters
OR-36790	IV POLE SUSPENSION POLE AUTO-CLUTCH adjustable height from 33" - 57" Recommended for ceiling up to 11 ft. w

HOOKS & RAM HORNS

HEAVY-DUTY RAM HORNS, U-HOOKS OR A COMBINATION OF BOTH AVAILABLE



COLOUR-CODING OFFERED FOR QUICK IDENTIFICATION BY DEPARTMENT

POLES & HANDLES

EASILY CLEANED AND DISINFECTED. RESISTS RUST AND PITTING



THE LARGE PRONG KNOB FOR HEIGHT ADJUSTMENT IS EASY TO TURN EVEN WITH WET HANDS

BASES & CASTERS

OPTIONAL WEIGHT KIT AVAILABLE FOR INCREASED STABILITY




SMOOTH-ROLLING NON-MARKING DURABLE CASTERS



HOOKS & RAM HORNS

ITEM #	DESCRIPTION
OR-36499	STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 2
OR-36500	STAINLESS STEEL RAM HORNS, AND ENDCAP Set of 4
OR-36501	STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 2
OR-36502	STAINLESS STEEL U-HOOKS, AND ENDCAP Set of 4
OR-36504	RAM HORNS (2), U-HOOKS (2), AND ENDCAP Set of 4
OR-36503H	RAM HORN SET/SUSPENSION POLE
OR-36503U	U-HOOK SET/SUSPENSION POLE
V-50002	URINARY DRAIN HOOK, AND COLLAR WITH KNOB Set of 1

*U-Hook & Ram Horn sold separately



OR-36504



OR-36499

HANDLES



V-50003 BLACK STEEL
ROUND HANDLE



V-50006 STAINLESS STEEL
TRIANGULAR HANDLE



V-50007 CREAM
SEMI-ROUND HANDLE

BASES & CASTERS

ITEM #	DESCRIPTION
V-12007	CASTERS 2" GRIP RING STEM OR-36770. Used on V-40016
V-12001	CASTERS 2"/HOLE/GREY OR-36760-25,OR-36771-25
V-12003	CASTERS 3" OR-36760-253,OR-36771-253
V-40013-25	STAINLESS STEEL BASE 25" base. Sold without casters
V-40016	ALUMINIUM BASE, CREAM 23" base. Sold without casters



V-40014-25-ST



V-40013-25



V-40016