

# INSTRUMENT/ SUPPLY CABINET

Our Instrument/Supply Cabinets are designed to provide enclosed storage for O.R. supplies, medications, instrumentation, and apparatus.

All Instrument/Supply Cabinets are backed by our **1 year warranty** on structural and all components.

## KEY FEATURES

- Cabinets come with four easily adjustable shelves (on 1" / 25 mm increments) and are fully interchangeable with other instrument/supply cabinets of similar width
- Standard sliding glass doors are framed in stainless steel and provide quick access to cabinet contents (optional door systems available)
- Cabinets are 22" / 559 mm deep x 84" / 2134 mm high at front
- Optional sloped top section adds 9 ½" / 241 mm to overall cabinet height
- Optional depth cabinets of 17 ¾" / 450 mm are available upon request



OR-3782

## PRODUCT SELECTION

ITEM #	DESCRIPTION
<b>OR-3782-A</b>	INSTRUMENT/SUPPLY CABINET Width: 47" (1194 mm)
<b>OR-3782-B</b>	INSTRUMENT/SUPPLY CABINET Width: 35" (889 mm)
<b>OR-3782-C</b>	INSTRUMENT/SUPPLY CABINET Width: 24" (610 mm)
<b>OR-3782-GH</b>	HINGED GLASS DOORS
<b>OR-3782-RF</b>	RECESSED FLANGE KIT
<b>OR-3782-ST</b>	SLOPED TOP SECTION