


# INSTRUMENT TABLES

Our Instrument Tables come in six versatile models and are an economical back table for instruments or dressings. Fully welded heavy gauge stainless steel construction polished to a #4 finish provides excellent rigidity and extended product life.

All instrument tables are backed by an outstanding **7 year warranty** on structural and all components (excluding casters).

## KEY FEATURES

- Channel support sub-frame for maximum weight capacity and structural rigidity
- Fully welded seamless tubular legs H-frame supports or lower shelf assembly
- Non-marking swivel design casters come in 2" (51 mm) or 3" (76 mm) depending on table size
- All models are 34" (864 mm) high and feature a lower shelf
- Open tubular base design is available upon request (reference # OR-3643)



## PRODUCT SELECTION

ITEM #	DIMENSIONS (L x W)
<b>OR-3641A</b>	20" x 16" (506 mm x 406 mm)
<b>OR-3641B</b>	30" x 16" (762 mm x 406 mm)
<b>OR-3641C</b>	34" x 18" (864 mm x 457 mm)
<b>OR-3641D</b>	36" x 20" (914 mm x 508 mm)
<b>OR-3641H*</b>	46" x 24" (1168 mm x 610 mm)
<b>OR-3641L*</b>	58" x 24" (1473 mm x 610 mm)

\*H and L models supplied with 3" / 76 mm caster assemblies

## ACCESSORIES

ITEM #	DESCRIPTION
<b>OR-3641X-R3</b>	THREE-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
<b>OR-3641X-R4</b>	FOUR-SIDED GUIDE RAIL Assembly for Table Top or Lower Shelf (Specify size)
<b>OR-3640</b>	DRAWER ASSEMBLY Inside Dimensions: 15" L x 16 3/4" W x 3 3/4" D (381 mm x 426 mm x 95 mm)