

KN95 DISPOSABLE PARTICULATE RESPIRATOR

The Batrik KN95 Disposable Particulate Respirator is designed to protect wearers from airborne particles. The **latex free** respirator is made with spunbond non-woven and melt blown layers and features an **ear-loop and adjustable nose bridge clip** designed to provide the ultimate in comfort when worn over a prolonged period.



Performance Standard: GB2626 - 2006

Filtration Efficiency: ≥ 95%

Certifications: Bacterial filtration efficiency ≥ 95%.Inhalation and exhalation resistance tested and conform to NIOSH criterion (NIOSH procedure TEB-APRSTP-0007 & TEB-APR-STP-0003)

FEATURES

- Secured seal: adjustable nose piece for optimize seal. Advanced non-woven fabric makes it easy to breathe.
- Ergonomic Design: durable & stretchable straps with dual point attachments, ensuring reliable protection.

PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
SF-KN95	KN95 Particulate Respirator	40 / Case