

MAYO STANDS

Our Mayo Stands come in seven models and are a great choice for general clinical use and operating room applications. The stainless steel seamless tubular construction is fully welded and polished for ease of cleaning and maximum structural support.

All mayo stands are backed by an outstanding **7 year warranty** on structural and all components (excluding casters).


KEY FEATURES

- Auto-clutch height adjustment system automatically locks tray at desired height and provides for smooth raising and lowering
- Safety mechanism prevents the tray frame from unexpectedly dropping
- Casters are ball bearing swivel 2" (51 mm) diameter
- All models feature a tray height adjustment range of 38" to 58" (965 mm to 1473 mm)

STAINLESS STEEL TRAYS ARE POLISHED, AND DIE-FORMED WITH RADIUS CORNERS. THEY COME IN TWO STANDARD SIZES

PRODUCT SELECTION

ITEM #	DESCRIPTION	TRAY SIZE (L X W)
OR-3657	EMERGENCY MAYO STAND WITH 4 CASTER BASE, AND T-LOCK Manually Operated	19"x13" (483 mm x 330 mm)
OR-3658	STANDARD TRAY Push-Button Sled-Type Base and Auto-Clutch System	19"x13" (483 mm x 330 mm)
OR-3662	LARGE TRAY Push-Button Sled-Type Base and Auto-Clutch System	21"x16" (533 mm x 406 mm)
OR-3660	STANDARD TRAY Manually Operated with Two Casters Sled-Type Base and T-Lock Handle	19"x13" (483 mm x 330 mm)
OR-3661	LARGE TRAY Manually Operated with Two Casters Sled-Type Base and T-Lock Handle	21"x16" (533 mm x 406 mm)
OR-3665	STANDARD TRAY Foot-Operated Mayo Stand Sled-Type Base and Foot Release Auto-Clutch System	19"x13" (483 mm x 330 mm)
OR-3666	LARGE TRAY Foot-Operated Mayo Stand Sled-Type Base and Foot Release Auto-Clutch System	21"x16" (533 mm x 406 mm)



OR-3665