


BATRIK PIN GUARDS



- Safely cover the sharp ends of your K-Wire and Steinmann Pins
- Color coded by wire size which eliminates the risk of selecting of the wrong size
- A removable sleeve is provided to avoid injury and facilitate application (sterile only)
- Non-sterile model is available with or without plastic sleeves
- Packaging designed to reduce waste
- Latex free



PRODUCT SELECTION

ITEM # STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS + SLEEVE	ITEM # NON-STERILE PIN GUARDS (NO SLEEVE)	COLOR	DIMENSIONS (WIRE SIZE)
14-1028	914-1028 - SLV	914-1028	Yellow	0.028" / 0.71 mm
14-1035	914-1035 - SLV	914-1035	Light Blue	0.035" / 0.89 mm
14-1045	914-1045 - SLV	914-1045	White	0.045" / 1.14 mm
14-1054	914-1054 - SLV	914-1054	Dark Red	0.054" / 1.40 mm
14-1062	914-1062 - SLV	914-1062	Green	0.062" / 1.57 mm
14-1078	914-1078 - SLV	914-1078	Pink	0.078" / 1.98 mm
14-1093	914-1093 - SLV	914-1093	Grey	0.093" / 2.36 mm
14-1109	914-1109 - SLV	914-1109	Beige	0.109" / 2.77 mm
14-1125	914-1125 - SLV	914-1125	Purple	0.125" / 3.17 mm
14-1141	914-1141- SLV	914-1141	Brown	0.141" / 3.58 mm
14-1156	914-1156 - SLV	914-1156	Dark Blue	0.156" / 3.96 mm
14-1188	914-1188 - SLV	914-1188	Red	0.188" / 4.77 mm

STERILE PACKAGING: 80 units/box (2 units per sterile pouch - 40 pouches per box)

NON-STERILE PACKAGING: 100 units/box