

Batrik PrAmLi Power³ LUBE™ Instrument Lubricant

Specially formulated to lubricate surgical instruments, hinges, box locks, and general equipment during autoclaving. PrAmLi Power³ LUBE™ is uniquely formulated with an all natural corrosion inhibitor.

- Forms a protective coating which can extend the instrument's, useful life
- Bacteriostatic
- Neutral pH 6.5 - 7.5
- Prevents corrosion

TIP:

FOR MILK BATH
SOLUTION, DILUTE
1 PART SOLUTION
IN 6 PARTS
WATER



Formulated with
PrAmLi Power³ Enzymes™



PRODUCT SELECTION

ITEM #	DILUTION LUBRICANT / WATER	PACKAGING	CASE DIMENSIONS
SS-21-6001	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	1 x 3.8 L (1 Gallon)	6" x 6" x 13"
SS-21-6004	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	4 x 3.8 L (4 Gallons)	13" x 13" x 14"
SS-21-6025	1 to 4 ml / L (0.125 to 0.5 oz / Gallon)	2 x 9.5 L (2 x 2.5 Gallons)	14.5" x 9.25" x 15"
SS-21-6032	No dilution required	1 x 946 ml (32 oz) Spray bottle	10" x 4"
SS-21-6012	No dilution required	12 X 946 ml (32 oz) Spray Bottles	18" x 12" x 11"