


STACKABLE STEP STOOLS

Made of lightweight cast aluminium and protected with an antimicrobial powder coating, our patented Stackable Step Stools have been designed with infection prevention as a prime focus.



PATENTED

- Weight of stool: 8.3 lbs (3.8 kg)
- Dimensions: 17.4" L x 13" W x 4.7" H (44.2 cm L x 33 cm W x 11.9 cm H)
- Antimicrobial coating prevents the spread of dangerous micro-organisms and suppresses the growth of mold, bacteria, yeast, and mildew
- Cart washable. Model OR-36364 (without mat) can be sterilized
- Anti-slip surface for enhanced stability
- Rubber-free corners avoids fluid collection
- No welds or large hollow openings to trap debris or fluids
- Transport cart available
- 300 lb (136 kg) weight capacity
- 7 year warranty covers all parts (step stool only)

PRODUCT SELECTION

ITEM #	DESCRIPTION	COLOR
OR-36363	Stackable Step Stool with mat	Black
OR-36363-02	Stackable Step Stool with mat	Dark Blue
OR-36363-03	Stackable Step Stool with mat	Green
OR-36363-07	Stackable Step Stool with mat	Orange
OR-36364	Stackable Step Stool without mat	N/A

Packaging: Sold by Unit

OPTIONAL ACCESSORIES

ITEM #	DESCRIPTION	PACKAGING
OR-3693	Transport cart	1 Unit
809-1705	Reusable Anti-Fatigue Mat with adhesive strips to secure it in place	5 / Box
899-1710	Single-Use Anti-Fatigue Mat with adhesive strips to secure it in place	10 / Box

REPLACEMENT PARTS

ITEM #	DESCRIPTION	PACKAGING
OR-36363-1A	Replacement Mat for OR-36363	1 Unit
OR-36363-1B	Replacement Foot Pegs	1 Unit

TIP: COLOR-CODE YOUR STEP STOOLS BY DEPARTMENT

