

STAINLESS STEEL MATERNITY BASSINET

Our stainless steel maternity bassinet is an economical bassinet for use in the operating room and/or maternity suite. The fully welded tubular frame construction provides maximum strength and facilitates easy cleaning and extended product life.

Our stainless steel maternity bassinet is backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

KEY FEATURES

- Modular frame allows for several storage designs including Trendelenburg positioning from either end of the frame
- Clear PVC infant basket features air/drainage holes and card holder
- Basket mattress is wrapped in a fire-retardant Med-Stat antibacterial covering
- Heavy duty 3" (76 mm) non-marking grey rubber casters (2 x Swivel and 2 x Total Lock)
- Lightweight design is easy to maneuver



TRENDELENBURG
POSITIONING
→



PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
PH-2655	TUBULAR FRAME MATERNITY BASSINET 33" x 17" x 38" (838 mm x 432 mm x 965 mm)
PH-2657	MATERNITY BASKET WITH CARD HOLDER High impact seamless acrylic clear Size: 27.5" x 13.5" x 9.5" / 698.5 mm x 342.9 mm x 241.3 mm Card Holder Size: 5" x 3.5"
PH-2659	BASKET MATTRESS Fire retardant anti-bacterial Size: 26" x 12.5" x 2" / 660.4 mm x 317.5 mm x 50.8 mm