

STRADDLE (NEUROSURGERY) TABLES

Our Straddle (Neurosurgery) Tables come in four lengths and are a convenient height-adjustable instrument table for cases involving many instruments and/or diagnostic equipment. Heavy-duty polished stainless steel construction with fully welded tubular support base.

All straddle (neurosurgery) tables are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- Direct-action drive mechanism provides smooth adjustment with no creeping or slippage
- Optional accessories include stainless steel extension shelves and marine-edge spill proof stainless steel work surface
- Feature four 3" (76 mm) grey non-marking casters (2 with swivel, 2 with brake)
- Height-adjustable from 36" to 58" (914 mm to 1473 mm) with stainless steel collapsible crank handle
- All models feature a standard 27.5" (700 mm) minimum clearance between the vertical leg supports



OR-3646X

Standard 27 1/2" / 700 mm minimum clearance between the vertical leg supports

PRODUCT SELECTION

| ITEM # | DESCRIPTION |
|-----------------|---|
| OR-3646X | STRADDLE TABLE, WITH STANDARD POSITION CRANK HANDLE Available lengths: A: 24" / 610 mm B: 32" / 813 mm C: 40" / 1016 mm D: 46" / 1168 mm |

Add suffix to Model # for work surface length required.
Standard Table / Extension Shelf width is 34" / 864 mm.

Stainless steel extension shelves and marine-edge work surfaces are available as optional features.

ACCESSORIES

| ITEM # | DESCRIPTION |
|--------------------|--|
| OR-3646ES | STAINLESS STEEL EXTENSION SHELF Available lengths: 4" and 8" / 102 mm and 203 mm depending on table width |
| OR-3646LH | LOWER POSITION CRANK HANDLE ASSEMBLY |
| OR-3646ME-X | MARINE-EDGE SPILL PROOF STAINLESS STEEL WORK SURFACE |