

UTILITY CARTS

Our Utility Carts are built to last using heavy-gauge grade 18-8 stainless steel and come in four versatile models for general purpose transport/storage use of supplies, instruments, or equipment in the Operating Room, trauma, emergency or minor procedure area.

All utility carts are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

SMART DESIGN

- Fully welded and polished to a #4 finish
- Spill-proof 1 inch (51 mm) lip on three sides of top surface and middle shelf
- Non-marking bumpers on push/pull handle and front corners of the cart
- **Reinforcing steel rod** is secured below the top shelf for **added strength and stability**
- **Steel plate** is affixed underneath all shelves to prevent wobble and bottoming-out

SUPERIOR MANEUVERABILITY

- Each cart features high load rating swivel casters that are mounted to heavy-duty chassis base frame
- Push/pull handle runs the full width of the cart and is positioned for ease of use



OR-3692F

PRODUCT SELECTION

ITEM #	SHELF SIZE (W x L)	SHELF CLEARANCE	OVERALL HEIGHT	CASTER DIAMETER	LOAD CAPACITY
OR-3692A	15 1/2" x 24" (394 mm x 610 mm)	12" (305 mm)	32" (813 mm)	3" (76 mm)	200 lbs (91 kg)
OR-3692B	17 3/4" x 30" (451 mm x 762 mm)	11 3/4" (298 mm)	32" (813 mm)	3" (76 mm)	200 lbs (91 kg)
OR-3692D	17 3/4" x 30" (451 mm x 762 mm)	12 1/8" (308 mm)	33 1/4" (845 mm)	4" (102 mm)	400 lbs (182 kg)
OR-3692F	20 3/4" x 35" (527 mm x 889 mm)	13 1/2" (343 mm)	37 1/8" (943 mm)	5" (127 mm)	500 lbs (224 kg)
OR-3692-WO*	19" x 31" (483 mm x 787 mm) Area 4.2 SQ FT (0.4 SQ M) Volume: 12 CU FT (0.3 CU M) Weight: 42 lbs. (19.1 kg)	Customizable	34" (864 mm)	4" (102 mm)	500 lbs (181 kg)

*Three sides have 1" high raised edge on top and middle shelf.

NOTE: Due to continual product development, specifications are subject to change without notice.