

# UTILITY INSPECTION TABLE

Our Utility Inspection Tables facilitate the assembly, inspection, and packing of instruments and/or equipment before or after the decontamination and sterilization process. All table components are built to last using heavy gauge Series 300 Grade 18-8 stainless steel and are die-formed, fully welded, and polished to a #4 blended finish.

All Utility Inspection Tables are backed by our **1 year warranty** on structural and all components.

## KEY FEATURES

- Rigid channel support design and tubular leg cross bracing ensures maximum stability, high load capacity, and extended product life
- Countertop is supported with die-formed stainless steel channels and integrally welded to the tubular leg frame base assembly
- Tables feature height adjustable stainless steel bullet feet
- Optional drawer assemblies, full length lower shelf, and overshelf assemblies provide flexible storage configurations
- Optional drawer assemblies are die-formed stainless steel with a double pan drawer front and heavy-duty ball bearing slide assemblies
- Custom sizes and designs are available upon request



SP-8724

## PRODUCT SELECTION

| ITEM #         | DESCRIPTION   |
|----------------|---|
| <b>SP-8710</b> | UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY<br>Length available: 44" / 1118 mm |
| <b>SP-8720</b> | UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY<br>Length available: 60" / 1524 mm |
| <b>SP-8730</b> | UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY<br>Length available: 76" / 1930 mm |
| <b>SP-8740</b> | UTILITY INSPECTION TABLE WITH TUBULAR LEG ASSEMBLY<br>Length available: 92" / 2337 mm |
| <b>SP-8712</b> | OVERSHELF ASSEMBLY WITH MOUNTING KIT<br>Length available: 44" / 1118 mm               |
| <b>SP-8722</b> | OVERSHELF ASSEMBLY WITH MOUNTING KIT<br>Length available: 60" / 1524 mm               |
| <b>SP-8732</b> | OVERSHELF ASSEMBLY WITH MOUNTING KIT<br>Length available: 76" / 1930 mm               |
| <b>SP-8742</b> | OVERSHELF ASSEMBLY WITH MOUNTING KIT<br>Length available: 92" / 2337 mm               |
| <b>SP-8714</b> | MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY<br>Length available: 44" / 1118 mm      |
| <b>SP-8724</b> | MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY<br>Length available: 60" / 1524 mm      |
| <b>SP-8734</b> | MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY<br>Length available: 76" / 1930 mm      |
| <b>SP-8744</b> | MODULAR DRAWER WITH FOUR SIDED APRON ASSEMBLY<br>Length available: 92" / 2337 mm      |