



MATERNITY BASSINETS



KERRI BEAR MATERNITY BASSINET™

The Kerri Bear Maternity Bassinet™ is easy to handle and clean, making it the ideal choice for hospitals looking for a versatile, safe, stylish, and long-lasting bassinet.

PRODUCT COMPONENTS

- Infant carrier
- Integrated changing table
- Spacious lower storage cabinet available in 5 colours: Designer White, Silverpine, Featherstone, Shadow Blue, and Mauve Blush

KEY FEATURES

ANTIMICROBIAL COATING

- Manufactured with our specialized lead-free, antibacterial protective powder paint coating*
- Prevents the spread of dangerous micro-organisms and suppresses the growth of mold and bacteria**
- Durable coating resists chipping and discoloration

BASSINET BASKET

- Frame and carrier are fully welded and painted with a baked protective anti-bacterial powder coating
- Carrier frame accommodates three basket positions: Neutral, Trendelenburg, and Reverse Trendelenburg
- Basket carrier can be easily adjusted back to reveal a convenient changing table underneath (pull safety switch)
- Safety switch mechanism locks basket carrier in place to prevent movement during transport

BASE CABINET

- Easy-pull handles and magnetic catches for easy door opening and secure door closing
- Fabricated from easy to disinfect high-quality materials (no wood or laminate components used)
- Large storage space for quick access to supplies and materials with optional shelf available



* Antibacterial paint specifications available upon request

** Laboratory documentation available upon request

MATTRESSES

- Two mattresses available: one for basket and one for changing table
- Mattresses made from fire retardant, antibacterial medium density foam material for infant support and comfort
- Changing table features a three-sided guard rail to ensure mattress is secured/prevent sliding



CASTERS

- Heavy duty 3" (76 mm) non-marking grey rubber casters facilitate excellent maneuverability of the bassinet
- 2 x Swivel and 2 x Total Lock








CASTERS

DIMENSIONS (L x W x H)

OVERALL SIZE	32.5" x 19.5" x 37.75" 826 mm x 495 mm x 959 mm	CHANGING TABLE MATTRESS	31" x 18.5" x 1" 787 mm x 470 mm x 25 mm
ACRYLIC BASKET	27.5" x 13.5" x 9.5" (Head End) 700 mm x 340 mm x 241 mm	BASKET MATTRESS	26" x 12.5" x 2" 660 mm x 318 mm x 50 mm
STORAGE CABINET (INTERIOR)	31" x 18" x 20" 787 mm x 457 mm x 508 mm		

PRODUCT SELECTION

ITEM #	COLOUR	
PH-2695WH	DESIGNER WHITE	
PH-2695SP	SILVERPINE	
PH-2695FS	FEATHERSTONE	
PH-2695SB	SHADOW BLUE	
PH-2695MB	MAUVE BLUSH	

ACCESSORIES

ITEM #	DESCRIPTION (L x W x H)
PH-2657	MATERNITY BASKET WITH CARD HOLDER High impact seamless acrylic clear Size: 27.5" x 13.5" x 9.5" / 698.5 mm x 342.9 mm x 241.3 mm Card Holder Size: 5" x 3.5
PH-2659	BASKET MATTRESS Fire retardant anti-bacterial Size: 26" x 12.5" x 2" / 660.4 mm x 317.5 mm x 50.8 mm
PH-2658	CHANGING TABLE MATTRESS Fire retardant anti-bacterial Size: 31" x 18.5" x 1" / 787.4 mm x 469.9 mm x 25.4 mm
PH-2672	INTERIOR CABINET WITH OPTIONAL SHELF Not adjustable
12008	CASTER 3" ALL SWIVEL All surface
12009	CASTER 3" WHEEL LOCK All surface, wheel lock
PH-2660	JAR RACK White only
PH-2661D	BASIN HOLDER, AND STAINLESS STEEL BASIN
PH-2671	CHART HOLDER

Find the images for the Items codes through the document



FOR EVERY NEW KERRI BEAR
BASSINET™ SOLD,
A DONATION OF \$10 WILL
BE MADE TO THE MISSING
CHILDREN'S NETWORK³.

3. Imperial Surgical has established special arrangements with the respective charity.

WE OFFER ONE OF THE BEST WARRANTIES IN THE INDUSTRY

All Imperial Surgical Kerri Bear Bassinets™ are warranted to be free from defects in workmanship and materials for a period of **seven years from date of original purchase**.

Accessories which are not manufactured by Imperial Surgical are warranted from defects in workmanship for one year.



STAINLESS STEEL MATERNITY BASSINET

Our stainless steel maternity bassinet is an economical bassinet for use in the operating room and/or maternity suite. The fully welded tubular frame construction provides maximum strength and facilitates easy cleaning and extended product life.

Our stainless steel maternity bassinet is backed by an outstanding **7 year warranty** on structural and all components (excluding casters) with one year labor included.

KEY FEATURES

- Modular frame allows for several storage designs including Trendelenburg positioning from either end of the frame
- Clear PVC infant basket features air/drainage holes and card holder
- Basket mattress is wrapped in a fire-retardant Med-Stat antibacterial covering
- Heavy duty 3" (76 mm) non-marking grey rubber casters (2 x Swivel and 2 x Total Lock)
- Lightweight design is easy to maneuver



PRODUCT SELECTION

ITEM #	DESCRIPTION (L x W x H)
PH-2655	TUBULAR FRAME MATERNITY BASSINET 33" x 17" x 38" (838 mm x 432 mm x 965 mm)
PH-2657	MATERNITY BASKET WITH CARD HOLDER High impact seamless acrylic clear Size: 27.5" x 13.5" x 9.5" / 698.5 mm x 342.9 mm x 241.3 mm Card Holder Size: 5" x 3.5
PH-2659	BASKET MATTRESS Fire retardant anti-bacterial Size: 26" x 12.5" x 2" / 660.4 mm x 317.5 mm x 50.8 mm



SURGMED GROUP

2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax: +1.514.631.9083
Email: info@surgmed.com | Website: www.imperialsurgical.com