



CADAVER STRETCHER

A #300 16 GAUGE STAINLESS STEEL
FABRICATION WITH FOUR 8" CASTERS AND
A HOLDING WEIGHT CAPACITY OF 300LBS.




CADAVER STRETCHER

Our Cadaver Stretcher is designed for careful and discreet transport. The stainless steel construction facilitates easy cleaning and decontamination and extended product life.

All cadaver stretchers are backed by our **1 year warranty** on structural and all components (excluding casters).

KEY FEATURES

- Features a lightweight removable canopy frame and vinyl covered canopy measuring 4" - 8" (10.2 cm – 20.3 cm)
- Complete circumference bumper provides protection
- 4 homa swiveling casters (2 x all swivel with directional lock, 2 x all swivel total break)



PRODUCT SELECTION

ITEM #	DESCRIPTION
PH-2753	CADAVER STRETCHER Stretcher and removable frame, vinyl covered canopy 8" (203 mm) Swivelling casters Holding weight capacity of 300 lbs
PH-2753HD	CADAVER STRETCHER - HEAVY DUTY Stretcher and removable frame, vinyl covered canopy 8" (203 mm) Swivelling casters Holding weight capacity of 650 lbs
12012	8" TOTAL BRAKES CASTERS With electrical resistance
12013	8" ALL SWIVEL DIRECTIONAL LOCK CASTERS With electrical resistance
13006-R0	COVER Thick resistant cover - Mint colour 75" L x 27" W x 17" D 1905 mm L x 686 mm W x 432 mm D

