



BRUSHES PORTFOLIO

Established in 1997. Proudly Canadian. 🇨🇦
Sold in over 80 countries worldwide.



BATRIK

ARC™ WIRE STYLE BRUSHES

SINGLE-USE AND REUSABLE

You've got options with our ARC™ Wire Brushes. Our single-use and reusable brushes are available in both flexible and rigid varieties and come in a multitude of sizes. ARC™ Wire Brushes are guaranteed to last longer than twisted wire and preserve the integrity of your valuable instruments.

BRUSH SELECTION TIPS:

INSTRUMENT SPECIFICATIONS

What is the maximum length of the instrument(s) to be cleaned?

Add at least 2" (5cm) to that length, so you have an adequate length to work with.

What is the inner diameter of the instrument?

Add at most 0.125" (3mm) to determine the diameter of the bristle to ensure effective cleaning. The bristle diameter has to slightly exceed the inner tube diameter. Bristles that

are too large will bend too much and not clean effectively.

BRUSH SELECTION

Simply find the section in our catalogue that best corresponds to the overall length and diameter of your instrument. As a rule, if your present brush is 8" (20.32cm) in length and we only offer 10" (25.4cm) length with your required diameter, that is entirely acceptable.


In many cases, our ARC™ Wire Brushes will have more bristles on the brush. As an example, many 5mm (0.19") diameter brushes are only 1" (2.5cm) in length. At BATRIK, we strived to provide the greatest amount of bristles to ensure the most effective cleaning.

CUSTOMIZABLE WIRE BRUSHES


A unique aspect of single-use ARC™ Wire Brushes is that they can be customized to any length:

- Up to 2000mm for Wire Brushes on plastic semi-rigid sheaths, rigid brushes from 5 to 20mm in diameter.
- No limit in length for Wire Brushes on flexible plastic sheaths, flexible brushes from 2-5mm in diameter

RIGID



FLEXIBLE



12" / 30 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-31502-S	41-31502	White	12.401" / 31.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-32503-S	41-32503	Yellow	12.795" / 32.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-32504-S	411-32504	Violet	12.795" / 32.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-32505-S	411-32505	Grey	12.795" / 32.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-35005-S	41-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35006-S	41-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-35008-S	41-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35010-S	41-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35012-S	41-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-35015-S	41-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-35020-S	41-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING Single-Use: Bag of 20 / Reusable: Bag of 3

18" / 45 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-46502-S	41-46502	White	18.307" / 46.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-47503-S	41-47503	Yellow	18.7" / 47.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-47504-S	411-47504	Violet	18.7" / 47.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-47505-S	411-47505	Grey	18.7" / 47.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-50005-S	41-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50006-S	41-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-50008-S	41-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50010-S	41-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50012-S	41-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-50015-S	41-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-50020-S	41-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single -Use: Bag of 20 / Reusable: Bag of 3

24" / 60 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-61502-S	41-61502	White	24.212" / 61.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1mm
	41-62503-S	41-62503	Yellow	24.6" / 62.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1mm
	411-62504-S	411-62504	Violet	24.6" / 62.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55mm
	411-62505-S	411-62505	Grey	24.6" / 62.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75mm
RIGID BRUSH & SHEATH	41-65005-S	41-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3mm
	41-65006-S	41-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3mm
	41-65008-S	41-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65010-S	41-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65012-S	41-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-65015-S	41-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-65020-S	41-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

30" / 75 CM SINGLE-END WIRE BRUSH

	ITEM # SINGLE-USE	ITEM # REUSABLE	COLOR	TOTAL LENGTH	BRUSH (Ø)	BRUSH LENGTH	SHEATH (Ø)
FLEXIBLE BRUSH & SHEATH	41-76502-S	41-76502	White	30.118" / 76.5 cm	0.078" / 2 mm	0.59" / 15 mm	0.39" / 1 mm
	41-77503-S	41-77503	Yellow	30.511" / 77.5 cm	0.118" / 3 mm	0.98" / 25 mm	0.39" / 1 mm
	411-77504-S	411-77504	Violet	30.511" / 77.5 cm	0.157" / 4 mm	0.98" / 25 mm	0.61" / 1.55 mm
	411-77505-S	411-77505	Grey	30.511" / 77.5 cm	0.19" / 5 mm	0.98" / 25 mm	0.68" / 1.75 mm
RIGID BRUSH & SHEATH	41-80005-S	41-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80006-S	41-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
	41-80008-S	41-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80010-S	41-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80012-S	41-80012	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
	41-80015-S	41-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
	41-80020-S	41-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Single-Use: Bag of 20 / Reusable: Bag of 3

FLEXIBLE 2MM




FLEXIBLE 3MM



FLEXIBLE 4MM



FLEXIBLE 5MM



FLEXIBLE

Flexible brush diameter (Ø) ranging from 2 mm to 5 mm

RIGID 20MM



RIGID 15MM



RIGID 12MM



RIGID 10MM



RIGID 8MM



RIGID 6MM



RIGID 5MM



RIGID

Rigid brush diameter (Ø) ranging from 5 mm to 20 mm

BRUSH SELECTION SYSTEM (Actual Size)

2mm
0.078"



3 mm
0.118"



4 mm
0.157"



5 mm
0.196"



6 mm
0.236"



8 mm
0.314"



10 mm
0.394"



12 mm
0.472"



15 mm
0.591"



20 mm
0.787"



BATRIK

REUSABLE DOUBLE ARC™ WIRE BRUSHES

ASYMMETRIC

14" / 35 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35056	White	13.779" / 35 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35810	Blue	13.779" / 35 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35112	Yellow	13.779" / 35 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35152	Green	13.779" / 35 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50056	White	19.685" / 50 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50810	Blue	19.685" / 50 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50112	Yellow	19.685" / 50 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50152	Green	19.685" / 50 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65056	White	25.59" / 65 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65810	Blue	25.59" / 65 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	25.59" / 65 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65152	Green	25.59" / 65 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM ASYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH 2 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80056	White	31.496" / 80 cm	0.19" / 5 mm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80810	Blue	31.496" / 80 cm	0.31" / 8 mm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-80112	Yellow	31.496" / 80 cm	0.39" / 10 mm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80152	Green	31.496" / 80 cm	0.59" / 15 mm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

BATRIK REUSABLE DOUBLE ARC™ WIRE BRUSHES SYMMETRIC

14" / 35 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-35005	White	13.779" / 35 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-35006	Green	13.779" / 35 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-35008	Blue	13.779" / 35 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-35010	Red	13.779" / 35 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-35012	Yellow	13.779" / 35 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-35015	Green	13.779" / 35 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-35020	Blue	13.779" / 35 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

20" / 50 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-50005	White	19.685" / 50 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-50006	Green	19.685" / 50 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-50008	Blue	19.685" / 50 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-50010	Red	19.685" / 50 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-50012	Yellow	19.685" / 50 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-50015	Green	19.685" / 50 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-50020	Blue	19.685" / 50 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

26" / 65 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-65005	White	25.59" / 65 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-65006	Green	25.59" / 65 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-65008	Blue	25.59" / 65 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-65010	Red	25.59" / 65 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65012	Yellow	25.59" / 65 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-65015	Green	25.59" / 65 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-65020	Blue	25.59" / 65 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

31" / 80 CM SYMMETRIC DOUBLE-END WIRE BRUSHES

ITEM #	COLOR	TOTAL LENGTH	BRUSH 1 (Ø)	BRUSH LENGTH	SHEATH (Ø)
40-80005	White	31.496" / 80 cm	0.19" / 5 mm	1.968" / 50 mm	0.118" / 3 mm
40-80006	Green	31.496" / 80 cm	0.23" / 6 mm	1.968" / 50 mm	0.118" / 3 mm
40-80008	Blue	31.496" / 80 cm	0.31" / 8 mm	1.968" / 50 mm	0.157" / 4 mm
40-80010	Red	31.496" / 80 cm	0.39" / 10 mm	1.968" / 50 mm	0.157" / 4 mm
40-65112	Yellow	31.496" / 80 cm	0.49" / 12.5 mm	1.968" / 50 mm	0.157" / 4 mm
40-80015	Green	31.496" / 80 cm	0.59" / 15 mm	1.968" / 50 mm	0.228" / 5.8 mm
40-80020	Blue	31.496" / 80 cm	0.78" / 20 mm	1.968" / 50 mm	0.228" / 5.8 mm

PACKAGING: Bag of 3

ASYMMETRIC BRUSHES

5 MM

6 MM



8 MM

10 MM



10 MM

12.5 MM



15 MM

20 MM




BATRIK

REUSABLE IN-BRUSH™ ANTIMICROBIAL TWISTED WIRE BRUSHES

Our In-Brush™ Antimicrobial Reusable Twisted Wire Brushes are designed for use with Frazier/Baron suction tubes. The twisted wires have a fine/medium diameter yet are rigid enough to ensure excellent cleaning capabilities.

Our wide assortment of Nylon Twisted Wire Brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

GOOD TO KNOW: OUR TWISTED
WIRE BRUSHES FEATURE
A PROTECTIVE BLUE TIP
FOR ADDED SAFETY.



10" / 250 MM FRAZIER /BARON SUCTION TUBE BRUSHES

ITEM #	SIZE	BRUSH LENGTH	BRUSH (Ø)
21-10039	3 FR	0.5" / 12.7 mm	0.039" / 1 mm
21-10067	5 FR	0.5" / 12.7 mm	0.065" / 1.67 mm
20-10059	5 FR	1" / 25.4 mm	0.059" / 1.50 mm
20-10079	6 FR	1" / 25.4 mm	0.079" / 2 mm
20-10090	7 FR	1" / 25.4 mm	0.090" / 2.3 mm
20-10106	8 FR	1" / 25.4 mm	0.106" / 2.7 mm
20-10118	9 FR	1" / 25.4 mm	0.118" / 3 mm
20-10129	10 FR	1.5" / 38 mm	0.129" / 3.3 mm
20-10146	11 FR	1.5" / 38 mm	0.0146" / 3.7 mm
20-10157	12 FR	1.5" / 38 mm	0.157" / 4 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

12" / 305 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
22-12078	2" / 50.8 mm	0.078" / 2 mm
22-12118	2" / 50.8 mm	0.118" / 3 mm
22-12157	3" / 76 mm	0.157" / 4 mm
22-12197	3" / 76 mm	0.197" / 5 mm
22-12236	3" / 76 mm	0.236" / 6 mm
22-12393	3" / 76 mm	0.394" / 10 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "- 50".

16" / 406 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
23-16197	2.5" / 63.5 mm	0.197" / 5 mm
23-16236	2.5" / 63.5 mm	0.236" / 6 mm
23-16276	3" / 76 mm	0.276" / 7 mm
23-16394	3" / 76 mm	0.394" / 10 mm
23-16472	3" / 76 mm	0.472" / 12 mm
23-16591	3" / 76 mm	0.591" / 15 mm
23-16787	3" / 76 mm	0.787" / 20 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

18" / 457 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
19-18010	0.47" / 12 mm	0.039" / 1 mm
19-18015	0.5" / 12 mm	0.06" / 1.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

24" / 610 MM NYLON TWISTED WIRE BRUSHES

ITEM #	BRUSH LENGTH	BRUSH (Ø)
24-24098	2.5" / 63.5 mm	0.098" / 2.5 mm
24-24197	2.5" / 63.5 mm	0.197" / 5 mm
24-24276	3" / 76 mm	0.276" / 7 mm
24-24394	3" / 76 mm	0.394" / 10 mm
24-24472	3" / 76 mm	0.472" / 12 mm


PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

TRACHEOSTOMY BRUSHES

ITEM #	DESCRIPTION	TOTAL LENGTH	BRUSH LENGTH	BRUSH (Ø)
25-08500	Large 6-12.3 3 / 4"	8" / 203 mm	3.75" / 95 mm	0.5" / 12.7 mm
25-37550	Small 000-5.1 7 / 8"	4" / 93 mm	1.88" / 48 mm	0.5" / 12.7 mm

PACKAGING: Bag of 3. Also available in bags of 50, add the suffix "-50".



BRONCHO BRUSHES

ITEM #	DESCRIPTION	BRUSH LENGTH	BRUSH (Ø) INCHES / MM
26-06175	Small	6" / 152 mm	1.75" / 44.5 mm
26-08325	Large	8" / 203 mm	3.25" / 82.6 mm

PACKAGING: Bag of 3.

Also available in bags of 50, add the suffix "-50".

BATRIK

STAINLESS STEEL ARC™ MICRO BRUSHES

Our Stainless Steel ARC™ Micro Brushes are some of the smallest brushes in the industry. They are designed to clean channels with an inner diameter ranging from 0.5 mm to 1.40 mm.

- Flexible stainless steel sheath will not bend like the regular metal twisted wire brushes
- Polyamide bristles prevents or slows bacterial growth between sterilizations
- Reusable*
- Available in two lengths

TIP: USING THE PULL-THROUGH
TECHNIQUE IS HIGHLY
RECOMMENDED WITH THESE
BRUSHES. SIMPLY INSERT THE
BLUNT, NON-BRUSH END FIRST.



PRODUCT SELECTION

ITEM #	TOTAL LENGTH	BRUSH SIZE (Ø)	BRUSH LENGTH	SHEATH (Ø)	PACKAGING
35-03015	12" / 30 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube
35-04515	18" / 45 cm	0.06" / 1.5 mm	0.6" / 15 mm	0.02" / 0.40 mm	6 / Tube

*Although all materials used in the fabrication of these brushes can be re-processed repeatedly, due to the delicate nature of these brushes re-use is conditional on their handling and application.

BATRIK

INSTRUMENT CLEANING FLAT BRUSHES

Our Instrument Cleaning Flat Brushes are specifically designed for hospital use in the pre-clean and preparation of instruments.

- Long ergonomic handle permits ease-of-use
- Extra long length of brush head (60 [±5] mm in Single-Use and 75 [±5] mm in Reusable)
- Transparent brush fibers allow for quick verification of its cleanliness
- The plastic handle will not release any toxic components
- Available in 6 models (4 reusable/autoclavable and 2 single-use)



FLEXIBLE

This brush is designed for large instruments that require a more gentle cleaning application such as optics, plastics, or hard to reach spot.



EXTRA FLEXIBLE

The very long, flexible bristles eliminate any stress during sanitation, and makes this the ideal brush for fragile instruments as well as both pliable and rigid endoscopes. The length of the bristles assist in cleaning areas that would otherwise be difficult to reach, such as the accessory channel or tip deflection of pliable endoscopes.



RIGID

The rigidity of this brush, along with the length of the bristles and handle, improves and simplifies instrument cleaning.



EXTRA RIGID

Instrument safe polyamide fibers offer steel-like rigidity and allows deep cleaning of soiled instruments, particularly in grooved or uneven surfaces.



50-40919



50-40921



50-40920



50-40918



50-40922



50-40917

PRODUCT SELECTION

ITEM #	BRUSH TYPE	DESCRIPTION	COLOR	BRUSH LENGTH	BRUSH SIZE	PACKAGING
50-40917	Extra Rigid	Reusable Autoclavable Flat Brush	Blue	220 mm	75 mm x 20 mm x 10 mm	Bag of 3
50-40918	Rigid	Reusable Autoclavable Flat Brush	Green	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40919	Flexible	Reusable Autoclavable Flat Brush	Yellow	220 mm	75 mm x 20 mm x 15 mm	Bag of 3
50-40920	Extra Flexible	Reusable Autoclavable Flat Brush	Red	220 mm	75 mm x 20 mm x 25 mm	Bag of 3
50-40921	Flexible	Single-Use Flat Brush	Orange	220 mm	60 mm x 20 mm x 15 mm	Bag of 20
50-40922	Rigid	Single-Use Flat Brush	Light Blue	220 mm	60 mm x 20 mm x 15 mm	Bag of 20


Reusable models are also available in bags of 50. Add the suffix "- 50".

BATRIK

REUSABLE ORTHOPEDIC RASPS CLEANING BRUSH

Our Reusable Orthopedic Rasps Cleaning Brush is used to clean orthopedic rasps and other similar instruments used in orthopedic surgery.

- Long lasting and can withstand many sterilization cycles
- Ergonomic steamable thermoplastic handle
- Antibacterial polyamide bristles
- Reusable
- Latex free



PRODUCT DESCRIPTION

ITEM #	TOTAL LENGTH	HANDLE LENGTH	BRUSH SIZE	PACKAGING
50-RASP	9" / 23 cm	5.5" / 14 cm	0.78" x 2.7" x 3.5" 2 cm x 7 cm x 9 cm	3 / Bag


BATRIK

REUSABLE PRECISION BRUSH

INSTRUMENT CLEANING

Our Reusable Precision Instrument Cleaning Brush is the ideal brush for the quick removal of difficult debris. It cleans externally as well as deep into instrument jaws. Although the firm and durable bristles feel like stainless steel, there is absolutely no risk of instrument damage or personal injury.

- Bristles have an extra firm density and superior durability
- Polyamide (TYNEX®) bristles
- Steamable thermoplastic ergonomic handle
- Steam / chemical sterilizable
- Latex free



GOOD TO KNOW: ELIMINATES
THE NEED FOR THE STAINLESS
STEEL TOOTHBRUSH-STYLE BRUSH,
WHICH CAN DAMAGE VALUABLE
INSTRUMENTS!

PRODUCT SELECTION

ITEM #	DESCRIPTION	TOTAL LENGTH
50-40800	Reusable Precision External Cleaning Brush	6.29" / 160 mm
PACKAGING: Bag of 3		




BATRIK

REUSABLE INSTRUMENT CLEANING BRUSHES

Designed specifically for hospital applications, our Reusable Instrument Cleaning Brushes feature long and flexible bristles to maximize the effectiveness of large-surface cleaning.

- High quality polyamide clear bristles
- Easy to identify debris
- Autoclavable
- Latex free

REUSABLE JAR/BOTTLE BRUSH**STERILIZER DRAIN BRUSH**

PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	BRUSH LENGTH	BRUSH (Ø)	TOTAL LENGTH
50-10295	Reusable Jar / Bottle Brush	120 mm ± 2 mm	120 mm ± 5 mm	7 mm ± 2 mm	13.19" / 355 mm
50-10224	Sterilizer Drain Cleaning Brush	N/A	4" / 10.1 cm	2" / 5.08 cm	24" / 609 mm
PACKAGING: Bag of 3					

BATRIK

ACETABULAR REAMER BRUSH

Our double-ended cleaning brush is manufactured exclusively for cleaning acetabular bone reamers. The “Pig Tail” rounded style circular head brush is designed to effectively clean inside the reamer, while the rigid straight brush on the reverse side is used for cleaning individual cutting holes. Special fusion of the brush heads to the hand piece along with its high quality antimicrobial medical grade bristles ensures a durable, reusable, and highly effective cleaning tool.



STRAIGHT BRUSH



ROUND BRUSH



PRODUCT SELECTION

ITEM #	DESCRIPTION	HANDLE LENGTH	STRAIGHT BRUSH	ROUND BRUSH HEAD
50-ACETAB	Acetabular Reamer Cleaning Brush	4.72\" / 120 mm Ø 0.47\" / 12 mm	1.18\" / 30 mm Ø 0.3125\" / 8 mm	2\" / 50 mm Ø 1.75\" / 45 mm

PACKAGING: Bag of 3

NOTE: Specifications subject to small variances

BATRIK**ALL PURPOSE REUSABLE BRUSHES****Plastic Handle Brushes**

We offer a wide selection of brushes for instrument cleaning applications. These brushes are ideal for large bore, uniquely shaped, and endoscopic instruments.

In-Brush™ Plastic Handle

Available in regular toothbrush style (nylon or stainless) and double-ended brushes.

**28-72550****28-69394****29-72550****29-69394****29-70050****50-40911****PRODUCT SELECTION**

ITEM #	DESCRIPTION	LENGTH TOTAL	BRUSH LENGTH	BRUSH HEIGHT
28-72550	Plastic Single End Nylon Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-72550	Plastic Single End Stainless Steel Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
29-70050	Plastic Single End Brass Reusable Handle Brush	7.25" / 184 mm	1.5" / 38 mm	0.5" / 12.7 mm
28-69394	Plastic Double End Nylon Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
29-69394	Plastic Double End Stainless Steel Reusable Handle Brush	6.9" / 175 mm	1.38" / 35 mm 1" / 25 mm	0.39" / 10 mm 0.2" / 5 mm
50-40911	Hand Scrub Nylon Brush Reusable autoclavable	4.09" / 104 mm	3.74" / 95 mm	1.1" / 28 mm

PACKAGING: Bag of 3 - For a bag of 50, add "-50" following the item code.

BATRIK

AUTOCLAVE CLEANING BRUSH KITS



ROUND BRUSH

- Features 3 holes for adjustable pole placement
- Adjustable brush placement (parallel or perpendicular) and longer brush length provides easy access to hard-to-reach corners
- Made with FDA approved polyester bristles

HAND BRUSH

- Ergonomic handle for better hand comfort
- Longer, FDA approved polyester bristles
- Easy to clean (no cavity)

PRODUCT SELECTION

REGULAR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVEREG	Autoclave REGULAR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 8.89 cm x 20.32 cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm

SUPERIOR AUTOCLAVE CLEANING BRUSH KIT

ITEM # (KIT)	DESCRIPTION	ITEM # (UNIT)	DESCRIPTION	DIMENSIONS
50-CLAVESUP	Autoclave SUPERIOR Brush Cleaning Kit	50-35900	1 Hand Brush	Length: 3.5" x 8" / 11.5 cm x 20.32cm
		50-55800	1 Rounded Corner Brush (Head Brush)	Head size: 5.5 x 10" / 14 x 25.4 cm
		50-EXTPOLE	1 Anodized Aluminum Extension Pole	2-4 ft / 60-122 cm
		SF-100	2 SAVEFACE™ Face Shields	9" / 23 cm

REPLACEMENT COMPONENTS



HAND BRUSH
50-35900

Packaging: 1



BRUSH HEAD
50-55800

Packaging: 1



EXTENSION POLE
50-EXTPOLE

Packaging: 1



SAVEFACE™
FACE SHIELDS SF-100

Packaging: 100 / Box

BATRIK

HYSTEROSONOGRAPHY CATHETER

Hystero-sonography (also known as sonohysterography) can now be performed safely, simply and atraumatically with the gentle infusion of saline into the uterus with the use of our H/S Catheter.


- Available in two sizes: 5FR and 7FR
- Minimally invasive and ideal for office use
- Gentle placement sheath included (decreasing the need for a tenaculum) to minimize patient discomfort
- Flexible and easy-to-use
- Total length of 410mm with a working length of 280 mm
- Latex Free

GOOD TO KNOW: ACHIEVE
IMPROVED IMAGERY
AND DIAGNOSIS WITH
SONOHYSTEROGRAMS

“We depend upon the reliability, precision and functionality of the Batrik H/S Catheter in our clinical and research operations.

Having performed more than 100,000 sonohysterography procedures, we can unequivocally say that this is the best HSG catheter that we have ever used


”
-Dr. Alex Hartman



PRODUCT SELECTION

ITEM #	DESCRIPTION	PACKAGING
51-5005	5 FR H/S Catheter	10 / box
51-5007	7 FR H/S Catheter	10 / box

100



SURGMED GROUP

2165 46th Ave., Lachine, Quebec H8T 2P1
Phone: +1 514-631-7988 | Toll Free: 800-661-5432 | Fax.: +1 514-631-9083
Email: info@surgmed.com | Website: www.batrik.com