

Your First ™

FIRST WIPES

Our cleanroom wipes are designed to do the first wipe-down of an endoscope after a procedure. They are made from hydroentangled non-woven material composed of 45% polyester and 55% cellulose. These wipes are manufactured and packaged in a cleanroom environment.



- Can be used for a cleanroom environment of ISO Class 5 through 8
- Production process guarantees extremely low particle generation and chemical extractables
- Compatible with virtually any solvent or solution for cleaning and wiping (such as **PrAmLi Power³ Enzymes™**)
- Double-bag packaging
- Available in two versatile sizes

PRODUCT SELECTION

ITEM #	DESCRIPTION	DIMENSIONS	PACKAGING	CASE DIMENSIONS
17-A0909	Non-woven Cleanroom Wipes	9" x 9" (23 cm x 23 cm)	300 / Bag 14 Bags / Case	18.5" x 10" x 25.25"
17-A1212	Non-woven Cleanroom Wipes	12" x 12" (31 cm x 31 cm)	150 / Bag 10 Bags / Case	13" x 13" x 19"