


BATRIK CORNER POCKETS™


With three models to choose from, our Corner Pockets™ portfolio is a truly comprehensive line of products designed to prevent tearing associated with sharp corners or feet on surgical trays.



18-CP02



18-CP07




18-CP06

| COLOUR | BLUE | ORANGE | PINK |
|-----------------------|---------------------------------|--|--|
| USAGE | SINGLE-USE | REUSABLE <i>withstands repeated sterilization up to 134°C</i> | REUSABLE <i>withstands repeated sterilization up to 134°C</i> |
| DIMENSIONS | 2.8" x 4" (7.3 cm x 10.3 cm) | 2.8" x 4" (7.3 cm x 10.3 cm) | 3.4" x 4.7" (8.5 cm x 12 cm) |
| LATEX-FREE | ✓ | ✓ | ✓ |
| NO MOISTURE RETENTION | ✓ | ✓ | ✓ |
| FABRICATION | Flexible, silicone-like plastic | Flexible, silicone-like plastic | Our most flexible plastic for easiest placement and most secure protection |

SOFT POCKETS™

Our Soft Pockets™ Foam Corner Protectors are made from high density foam to effectively cushion and protect the corners of surgical trays to prevent tearing. They are especially useful with heavier instrument sets.



PRODUCT SELECTION

| ITEM # | DESCRIPTION | PACKAGING |
|------------|---|---------------------------------|
| 18-CP02** | Corner Pockets™ <i>Single-Use, Blue</i> | 240 / Bag |
| 18-EVAL2** | Corner Pockets™ <i>Single-Use, Blue</i> | 24 / Bag |
| 18-CP07 | Corner Pockets™ <i>Reusable, Orange</i> | 80 / Bag |
| 18-CORNER | Corner Pockets™ <i>Reusable, Orange</i> | 320 / Box (4 Bags of 80 pcs) |
| 18-EVAL7 | Corner Pockets™ <i>Reusable, Orange</i> | 24 / Bag |
| 18-CP06 | Corner Pockets™ <i>Reusable, Pink</i> | 80 / Bag |
| 18-CPF-02 | Soft Pockets™ <i>Single-Use Foam Corner Protectors</i> | 80 / Box (24 Boxes / Pallet) |

** Although the single use Corner Pockets™ appear similar to the reusable models they are blended differently and have only been validated as single-use. Re-using them may compromise sterility.