



DISCOVER BATRIK

SHOOZ_{COVERS}™



AND PUT YOUR
BEST FOOT FORWARD

PROUDLY MADE IN CANADA



NEW!

SHOOZCOVERS™

Batrik ShoozCovers™ are manufactured in Canada from high quality 35 g/m² polypropylene fabric and 100% latex free elastics for performance and durability you can rely on.


Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

WHY YOU'LL LOVE THEM:

- **Low linting non-woven 35 g/m² polypropylene fabric** protects footwear and keeps floors clean
- Particle count tested and complies with **IEST-RP-CC003.4**
- Shooz-Repellent™ technology delivers **basic incidental fluid contact repellency**¹
- **Triple ultrasonic welds** increase tear strength and reduce linting
- High quality **made in the USA latex-free elastics** are triple reinforced to deliver a strong and secure fit
- **Generous sizing** and smart design for better fit and easy on-and-off convenience
- **Anti-skid tread on soles**
- **Truly 100% latex and silicone free shoe covers**²
- **Non-Sterile, Single Use**
- **Long-lasting to reduce waste**



MADE IN CANADA 



PERFECT FOR:

- Protecting a variety of footwear from shoes to clogs to boots
- Use in medical establishments, pharmaceutical labs, electronic manufacturing, GMP sites, and anywhere optimum hygiene is required
- Transitioning from “unclean” to “clean” areas

¹Shooz-Repellent™ technology provides basic incidental fluid contact repellency. It is not waterproof nor water resistant.

²Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

TWO VERSATILE SIZES:

UNIVERSAL


UN

Up to shoe size 11 (44.5 EU)
~ 42 cm x 16.5 cm / 16.5" x 6.5"³

EXTRA LARGE

XL

Shoe sizes 11-15 (44.5-49.5 EU)
~ 47 cm x 19 cm / 18.5" x 7.5"³



³SHOOZCovers™ are made longer (42 & 47cm) for accurate shoe size selection.

DONNING IS AS EASY AS 1-2-3:



MADE IN CANADA

- ☒ Support Canadian manufacturing
- ☒ Eliminate supply chain issues
- ☒ Reduce shipping costs

PRODUCT SELECTION

ITEM #	DESCRIPTION	SIZE	PACKAGING
SC-1216-PDB	ShoozCovers™ Disposable Shoe Covers <i>Blue 35 g/m² fabric</i>	Universal Up to shoe size 11 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	100 / Dispenser Box 6 Dispenser Boxes / Case
SC-1226-PDB	ShoozCovers™ Disposable Shoe Covers <i>Blue 35 g/m² fabric</i>	XL Shoe size 11 - 15 ~ 47 cm x 19 cm / 18.5" x 7.5"	100 / Dispenser Box 6 Dispenser Boxes / Case