



DISCOVER BATRIK

SHOOZCOVERS™

AND PUT YOUR BEST FOOT FORWARD



DISPOSABLE SHOE COVERS



PROUDLY MADE IN CANADA



SHOOZCOVERS™ | SOFT LINE™


Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality 35 g/m² polypropylene fabric and 100% latex free elastics for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.

ANTI-SKID TREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.



MADE IN CANADA 

OUR SOFT LINE™ IS PERFECT FOR:



Medical
Establishments



Pharmaceutical
Labs



GMP Sites



Electronic
Manufacturing

LOLLY BLUE | **SOFT LINE™**




WHY YOU'LL LOVE THE SOFT LINE™:

- Both shoe covers are manufactured in Canada with top-quality 35 g/m² polypropylene fabric and 100% latex-free USA-made elastics¹.
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request

UN XL

AVAILABLE IN TWO VERSATILE SIZES



BUBBLE GUM PINK | **SOFT LINE™**



SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!

UN

AVAILABLE IN UNIVERSAL



PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-1210-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)
SC-1216-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)
SC-1226-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 bags of 100 units)
SC-1220-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1000 (10 bags of 100 units)
SC-1616-PB	ShoozCovers™ - Soft Line™ 35 g/m ² polypropylene non-woven fabric	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)

ALTERNATIVE PACKAGING: For Items SC-1210-PB, SC-1220-PB add "-PDB" to end of code for a dispenser box instead of bags

¹Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

SHOOZCOVERS™ | TOUGH LINE™

EXTRA-LARGE FOR EXTRA-TOUGH WORK:


Our Batrik ShoozCovers™ Tough Line™ comes in an extra-large size that fits over a men's size 15 steel-toed boot, delivering optimal full coverage and protection while you work.¹

ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of ~55 g/m² polypropylene fabric, these highly durable covers provide added protection, increase tear strength and reduce linting.



Shoe sizes 12-15 (45-49.5 EU)
~ 47 cm x 19 cm / 18.5" x 7.5"



MADE IN CANADA 



OUR TOUGH LINE™ IS PERFECT FOR:



Homes or Commercial Buildings



Shielding against mud & grease



Medical Establishments



Pharmaceutical Labs

¹ ShoozCovers™ are made longer (42 & 47cm) for accurate shoe size selection.

BLACK LIQUORICE | TOUGH LINE™

WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit²
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request



XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:



Residential



Commercial



Industrial

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-1922-PB	ShoozCovers™ Tough Line™ ~55 g/m ² polypropylene non-woven fabric	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.

SHOOZCOVERS™ | COATED LINE™


STAY DRY:

Our Batrik ShoozCovers™ Coated Line™ is made from a ~55 g/m² blend of **polypropylene and polyethylene**, providing **high-quality water resistance**.¹ Ideal for environments where **maximum protection against fluids** is required.



TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with **triple ultrasonic welds**, a technique that joins materials together without the need for adhesives or stitching. This technique **increase tear strength and reduce linting**, for long lasting, resilient covers.



MADE IN CANADA 

OUR COATED LINE™ IS PERFECT FOR:



Medical Establishments



Pharmaceutical Labs



GMP Sites



Electronic Manufacturing

¹Our ShoozCovers™ Coated Line™ is resistant against all non petroleum based liquids.

MARSHMALLOW WHITE | COATED LINE™

WHY YOU'LL LOVE THE COATED LINE™:



- **Superior protection using our unique blend of fabrics** which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- **Elastics around the ankle create a strong and secure fit** which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- **Smart design** with a large opening, providing a better fit and easy on-and-off convenience
- **Other colors available upon request**



UN XL

AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

ITEM #	DESCRIPTION	COLOR	SIZE	PACKAGING
SC-2112-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 bags of 150 units)
SC-2112-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.88"	Case of 300 (2 boxes of 150 units)
SC-2122-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 bags of 120 units)
SC-2122-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 6.88"	Case of 240 (2 boxes of 120 units)

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags



2165 46th Ave., Lachine, QC, Canada, H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083
Email: info@surgmed.com | Website: www.batrik.com

Ver. 2023-10-26