



Disposable Shoe Covers
Soft Line™

DISCOVER BATRIK
SHOOZCOVERS™



Disposable Shoe Covers
Coated Line™



Disposable Shoe Covers
Tough Line™



Grippez Shooz™
Disposable Shoe Covers
Rhino Grey



Grippez Shooz™
Disposable Shoe Covers
Gecko Navy

ABOUT US

The **SURGMED Group of Companies** initially began its operations in 1992 as a medical device distribution company known as Advanced Surgi-Pharm Inc. Quickly becoming recognized as an industry leader serving the OBGYN, L&D, NICU, and Operating Room.

SURGMED made strategic acquisitions and greatly expanded into the direct manufacturing of medical devices under the brands of **BATRIK** and **IMPERIAL SURGICAL**. Today The SURGMED Group's North American manufactured medical devices are being sold in over 85 countries.


2165 46th Ave., Lachine, QC, Canada, H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083
Email: info@surgmed.com | Website: www.batrik.com



SHOOZCOVERS™ | SOFT LINE™

Batrik ShoozCovers™ Soft Line™ is manufactured in Canada from high quality **35 g/m² polypropylene fabric and 100% latex free elastics** for performance and durability you can rely on. Our anti-skid tread reduces the risk of slipping and our smart design delivers easy on-and-off convenience and great fit.

MADE IN CANADA



SHOOZ-REPELLENT™ TECHNOLOGY:

The **Soft Line™** integrates Shooz-Repellent™ technology which creates a protective barrier that repels fluids, delivering basic incidental fluid contact repellent. Shoe covers are not waterproof or water-resistant.



ANTI-SKID THREAD:

Our anti-skid tread reduces the risk of slipping with **triple ultrasonic welds** that increases tear strength and delivers easy on-and-off convenience.

OUR SOFT LINE™ IS PERFECT FOR:


MEDICAL
ESTABLISHMENTS

PHARMACEUTICAL
LABS


GMP SITES


ELECTRONIC
MANUFACTURING

LOLLY BLUE | **SOFT** LINE™



BUBBLE GUM PINK | **SOFT** LINE™

WHY YOU'LL LOVE THE SOFT LINE™:

- Both shoe covers are manufactured in Canada with top-quality 35 g/m² polypropylene fabric and 100% latex-free USA-made elastics¹.
- This line is particle count tested and complies with IEST-RP-CC003.4.
- Other colors available upon request.

SUPPORT A CAUSE, WEAR PINK:

- A portion of our Bubble Gum Pink shoe cover sales will be donated to charity to help further breast cancer research. Join us in making a meaningful contribution by choosing Pink ShoozCovers™!



PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-1210-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 1,000 (10 bags of 100 units)	20" x 16" x 16"
SC-1216-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25" x 8.875" x 18.75"
SC-1216-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5" x 19.5" x 14.25"
SC-1226-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i> <i>*Covers are available with or without skid-resistant tread.</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 600 (6 Dispenser Boxes of 100 units)	17.25" x 9.375" x 21.75"
SC-1220-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Lolly Blue	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 1,000 (10 bags of 100 units)	22" x 16" x 16"
SC-1616-PDB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 Dispenser Boxes of 100 units)	16.25" x 8.875" x 18.75"
SC-1616-PB	ShoozCovers™ - Soft Line™ <i>polypropylene non-woven fabric</i>	Bubble Gum Pink	Universal Up to shoe size 12 ~ 42 cm x 16.5 cm / 16.5" x 6.5"	Case of 600 (6 bags of 100 units)	19.5" x 19.5" x 14.25"

¹Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product.

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Softline	40	40 x 48 x 48 in	420 lbs	U	6 Dispenser Box / Case (100 un. each)

*Quantity of Cases per Pallet may vary

SHOOZCOVERS™ | TOUGH LINE™



MADE IN CANADA



EXTRA-LARGE FOR EXTRA-TOUGH WORK:

SHOE SIZES 12-15 (45-49.5 EU)
~ 47 CM X 19 CM / 18.5" X 7.5"

Our **Batrik ShoozCovers™ Tough Line™** comes in an extra-large size that fits over a **men's size 15 steel-toed boot**, delivering optimal full coverage and protection while you work.¹



ENHANCED DURABILITY:

These shoe covers are formulated with our **thickest fabric yet!** Made of **~55 g/m² polypropylene fabric**, these highly durable covers provide added protection, increase tear strength and reduce linting.

OUR TOUGH LINE™ IS PERFECT FOR:



HOMES OR
COMMERCIAL
BUILDINGS



SHIELDING AGAINST
MUD & GREASE



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS

¹ SHOOZCOVERSTM ARE MADE LONGER (42 & 47CM) FOR ACCURATE SHOE SIZE SELECTION.

BLACK LIQUORICE | TOUGH LINE™

WHY YOU'LL LOVE THE TOUGH LINE™:

- Triple ultrasonic welds that are long lasting to reduce waste
- High quality made in the USA latex-free elastics are triple reinforced to contract the shoe cover around your foot and deliver a strong and secure fit²
- Shooz-Repellent™ technology delivers basic incidental fluid contact repellency³
- Other colors available upon request


XL

AVAILABLE IN EXTRA LARGE

ANY TIME, ANYWHERE:


RESIDENTIAL



COMMERCIAL



INDUSTRIAL

PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-1922-PB	ShoozCovers™ - Tough Line™ Polypropylene non-woven fabric	Black Liquorice	XL Shoe size 12-15 ~ 47 cm x 19 cm / 18.5" x 7.5"	Case of 200 (2 bags of 100 units)	16" x 13.375" x 10"

ALTERNATIVE PACKAGING: Add "-PDB" to end of code for a dispenser box instead of bags.

PALLET INFORMATION:

	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Toughline	72	40 x 48 x 90 in	481 lbs	XL	2 Dispenser Box / Case (120 un. each)

*Quantity of Cases per Pallet may vary

² Many suppliers claim their products are latex free, however this often only relates to the fabric and not the elastics used. We source our ShoozCovers™ elastics from the USA to ensure that our components conform fully as a 100% latex free product. ³ Shooz-Repellent™ technology provides basic incidental fluid contact repellent. It is not waterproof nor water resistant.


SHOOZCOVERS™ | COATED LINE™

STAY DRY:

Our **Batrik ShoozCovers™ Coated Line™** is made from a **~55 g/m² blend of polypropylene and polyethylene**, providing high-quality water resistance.¹ Ideal for environments where maximum protection against fluids is required.



MADE IN CANADA



TRIPLE WELDED TECHNOLOGY:

These shoe covers are reinforced with triple ultrasonic welds, a technique that joins materials together without the need for adhesives or stitching. This technique increase tear strength and reduce linting, for long lasting, resilient covers.

OUR COATED LINE™ IS PERFECT FOR:



MEDICAL
ESTABLISHMENTS



PHARMACEUTICAL
LABS



GMP SITES



ELECTRONIC
MANUFACTURING

¹OUR SHOOZCOVERS™ COATED LINE™ IS RESISTANT AGAINST ALL NON PETROLEUM BASED LIQUIDS.

MARSHMALLOW WHITE | COATEDLINE™



WHY YOU'LL LOVE THE COATED Line™:

- Superior protection using our unique blend of fabrics which protects against not only water but all non-petroleum-based liquids, making them the ideal choice for environments where spills and leaks are a concern
- Elastics around the ankle create a strong and secure fit which is especially important when working in hazardous environments where a misplaced shoe cover could lead to contamination or injury
- Smart design with a large opening, providing a better fit and easy on-and-off convenience
- Other colors available upon request



UN - XL

AVAILABLE IN TWO VERSATILE SIZES

PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-2112-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-2112-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	150 / Dispenser Box 2 Dispenser Boxes / Case	16.62" x 14.875" x 8.75"
SC-2122-NPB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
SC-2122-NPDB	ShoozCovers™ - Coated Line™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Marshmallow White	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	120 / Dispenser Box 2 Dispenser Boxes / Case	16.62" x 14.875" x 8.75"

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
Coated	72	40 x 48 x 90 in	481 lbs	XL	2 Dispenser Box / Case (120 un. each)

*Quantity of Cases per Pallet may vary



GrippezShooz™

EXTRA-LARGE FOR EXTRA-TOUGH WORK:

GrippezShooz™ provide exceptional protection and **enhanced traction** for professionals working in the most demanding, slippery environments. Featuring waterproof protection, our shoe covers ensure safety and hygiene without compromising on quality. They are the ideal solution for tough jobs.



MADE IN CANADA



PERFECT FOR:



AEROSPACE



FOOD PROCESSING
INDUSTRY



OIL, GAS AND
MINING SITES



CONSTRUCTION AND
RENOVATION SITES

¹ INDUSTRIAL STRENGTH ELASTICS

GrippeezShooz™ ARE A PERFECT FIT FOR EVERYONE:

- Fast, Reliable Fit: Stretch material ensures a secure, easy-on fit. Coating provides protection against slippery surfaces.
- Comes in an extra-large size that fits over a men's size 15 steel-toed boot.
- Enhanced Safety: Anti-slip rubberized coating provides reliable grip on all surfaces, helping to prevent sliding in critical areas.
- Available in two colors: Gecko Navy and Rhino Grey.

RHINO GREY



GECKO NAVY



PRODUCT SELECTION:

ITEM#	DESCRIPTION	COLOR	SIZE	PACKAGING	CASE DIMENSIONS
SC-29312-NPDB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-29312-NPB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Rhino Grey	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-29322-NPDB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-29322-NPB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Rhino Grey	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"
SC-22212-NPDB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 dispenser boxes of 150 units)	16.62" x 14.875" x 8.75"
SC-22212-NPB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Gecko Navy	Universal Up to shoe size 12 ~ 42 cm x 17.5 cm / 16.5" x 6.9"	Case of 300 (2 bags of 150 units)	16" x 13.375" x 10"
SC-22222-NPDB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 dispenser boxes of 120 units)	16.62" x 14.875" x 8.75"
SC-22222-NPB	GrippeezShooz™ <i>Polypropylene + polyethylene non-woven fabric, waterproof</i>	Gecko Navy	XL Shoe size 12-15 ~ 47 cm x 17.5 cm / 18.5" x 6.9"	Case of 240 (2 bags of 120 units)	16" x 13.375" x 10"

PALLET INFORMATION:

PRODUCT	CASES PER PALLET*	PALLET DIMENSIONS	PALLET WEIGHT	SIZE	PACKAGING UOM
GRIPPEEZSHOOZ	50	40 x 48 x 88 in	446-479 lbs	U	2 Dispenser Box / Case (150 un. each)

*Quantity of Cases per Pallet may vary



2165 46th Ave., Lachine, QC, Canada, H8T 2P1
Phone: +1.514.631.7988 | Toll Free: 800.661.5432 | Fax.: +1.514.631.9083
Email: info@surgmed.com | Website: www.batrik.com